

DRAINAGE CATHETER SET




DRAINAGE CATHETER SET-DIRECT ACCESS



- Used for nephrostomy, abscess, and other multipurpose drainage.
- Catheter placement is performed using one step (direct access) technique.
- The safety string lock system provides strong fixation of the pigtail.
- Hydrophilic decreased friction of the distal catheter facilitates easy insertion.
- Large oval drainage holes increase drainage capability.
- Supplied sterile in peel-open packages.
- Intended for one-time use.



PIGTAIL TYPE (BT-PD1-series)



With SAFETY STRING LOCK						
Order Number	Size	Accepts Guide Wire Dia.	Working Lengh	Order Number	Size	Accepts Guide Wire Dia.
BT-PD1-0620-W(+CS)	6F	.035 inch (0.89mm)	19 cm (XX=20)	BT-PD1-1020-W(+CS)	10F	.038 inch (0.97mm)
BT-PD1-0625-W(+CS)	6F	.035 inch (0.89mm)	24 cm (XX=25)	BT-PD1-1025-W(+CS)	10F	.038 inch (0.97mm)
BT-PD1-0630-W(+CS)	6F	.035 inch (0.89mm)	29 cm (XX=30)	BT-PD1-1030-W(+CS)	10F	.038 inch (0.97mm)
BT-PD1-0635-W(+CS)	6F	.035 inch (0.89mm)	34 cm (XX=35)	BT-PD1-1035-W(+CS)	10F	.038 inch (0.97mm)
BT-PD1-0640-W(+CS)	6F	.035 inch (0.89mm)	39 cm (XX=40)	BT-PD1-1040-W(+CS)	10F	.038 inch (0.97mm)
BT-PD1-0645-W(+CS)	6F	.035 inch (0.89mm)	44 cm (XX=45)	BT-PD1-1045-W(+CS)	10F	.038 inch (0.97mm)
BT-PD1-0650-W(+CS)	6F	.035 inch (0.89mm)	49 cm (XX=50)	BT-PD1-1050-W(+CS)	10F	.038 inch (0.97mm)
BT-PD1-0720-W(+CS)	7F	.035 inch (0.89mm)	19 cm (XX=20)	BT-PD1-1220-W(+CS)	12F	.038 inch (0.97mm)
BT-PD1-0725-W(+CS)	7F	.035 inch (0.89mm)	24 cm (XX=25)	BT-PD1-1225-W(+CS)	12F	.038 inch (0.97mm)
BT-PD1-0730-W(+CS)	7F	.035 inch (0.89mm)	29 cm (XX=30)	BT-PD1-1230-W(+CS)	12F	.038 inch (0.97mm)
BT-PD1-0735-W(+CS)	7F	.035 inch (0.89mm)	34 cm (XX=35)	BT-PD1-1235-W(+CS)	12F	.038 inch (0.97mm)
BT-PD1-0740-W(+CS)	7F	.035 inch (0.89mm)	39 cm (XX=40)	BT-PD1-1240-W(+CS)	12F	.038 inch (0.97mm)
BT-PD1-0745-W(+CS)	7F	.035 inch (0.89mm)	44 cm (XX=45)	BT-PD1-1245-W(+CS)	12F	.038 inch (0.97mm)
BT-PD1-0750-W(+CS)	7F	.035 inch (0.89mm)	49 cm (XX=50)	BT-PD1-1250-W(+CS)	12F	.038 inch (0.97mm)
BT-PD1-0820-W(+CS)	8F	.038 inch (0.97mm)	19 cm (XX=20)	BT-PD1-1420-W(+CS)	14F	.038 inch (0.97mm)
BT-PD1-0825-W(+CS)	8F	.038 inch (0.97mm)	24 cm (XX=25)	BT-PD1-1425-W(+CS)	14F	.038 inch (0.97mm)
BT-PD1-0830-W(+CS)	8F	.038 inch (0.97mm)	29 cm (XX=30)	BT-PD1-1430-W(+CS)	14F	.038 inch (0.97mm)
BT-PD1-0835-W(+CS)	8F	.038 inch (0.97mm)	34 cm (XX=35)	BT-PD1-1435-W(+CS)	14F	.038 inch (0.97mm)
BT-PD1-0840-W(+CS)	8F	.038 inch (0.97mm)	39 cm (XX=40)	BT-PD1-1440-W(+CS)	14F	.038 inch (0.97mm)
BT-PD1-0845-W(+CS)	8F	.038 inch (0.97mm)	44 cm (XX=45)	BT-PD1-1445-W(+CS)	14F	.038 inch (0.97mm)
BT-PD1-0850-W(+CS)	8F	.038 inch (0.97mm)	49 cm (XX=50)	BT-PD1-1450-W(+CS)	14F	.038 inch (0.97mm)
BT-PD1-0920-W(+CS)	9F	.038 inch (0.97mm)	19 cm (XX=20)	BT-PD1-1620-W(+CS)	16F	.038 inch (0.97mm)
BT-PD1-0925-W(+CS)	9F	.038 inch (0.97mm)	24 cm (XX=25)	BT-PD1-1625-W(+CS)	16F	.038 inch (0.97mm)
BT-PD1-0930-W(+CS)	9F	.038 inch (0.97mm)	29 cm (XX=30)	BT-PD1-1630-W(+CS)	16F	.038 inch (0.97mm)
BT-PD1-0935-W(+CS)	9F	.038 inch (0.97mm)	34 cm (XX=35)	BT-PD1-1635-W(+CS)	16F	.038 inch (0.97mm)
BT-PD1-0940-W(+CS)	9F	.038 inch (0.97mm)	39 cm (XX=40)	BT-PD1-1640-W(+CS)	16F	.038 inch (0.97mm)
BT-PD1-0945-W(+CS)	9F	.038 inch (0.97mm)	44 cm (XX=45)	BT-PD1-1645-W(+CS)	16F	.038 inch (0.97mm)
BT-PD1-0950-W(+CS)	9F	.038 inch (0.97mm)	49 cm (XX=50)	BT-PD1-1650-W(+CS)	16F	.038 inch (0.97mm)

Without SAFETY STRING LOCK						
Order Number	Size	Accepts Guide Wire Dia.	Working Lengh	Order Number	Size	Accepts Guide Wire Dia.
BT-PD1-0520(+CS)	5F	.035 inch (0.89mm)	19 cm (XX=20)	BT-PD1-1020(+CS)	10F	.038 inch (0.97mm)
BT-PD1-0525(+CS)	5F	.035 inch (0.89mm)	24 cm (XX=25)	BT-PD1-1025(+CS)	10F	.038 inch (0.97mm)
BT-PD1-0530(+CS)	5F	.035 inch (0.89mm)	29 cm (XX=30)	BT-PD1-1030(+CS)	10F	.038 inch (0.97mm)
BT-PD1-0535(+CS)	5F	.035 inch (0.89mm)	34 cm (XX=35)	BT-PD1-1035(+CS)	10F	.038 inch (0.97mm)
BT-PD1-0540(+CS)	5F	.035 inch (0.89mm)	39 cm (XX=40)	BT-PD1-1040(+CS)	10F	.038 inch (0.97mm)
BT-PD1-0545(+CS)	5F	.035 inch (0.89mm)	44 cm (XX=45)	BT-PD1-1045(+CS)	10F	.038 inch (0.97mm)
BT-PD1-0550(+CS)	5F	.035 inch (0.89mm)	49 cm (XX=50)	BT-PD1-1050(+CS)	10F	.038 inch (0.97mm)
BT-PD1-0620(+CS)	6F	.035 inch (0.89mm)	19 cm (XX=20)	BT-PD1-1220(+CS)	12F	.038 inch (0.97mm)
BT-PD1-0625(+CS)	6F	.035 inch (0.89mm)	24 cm (XX=25)	BT-PD1-1225(+CS)	12F	.038 inch (0.97mm)
BT-PD1-0630(+CS)	6F	.035 inch (0.89mm)	29 cm (XX=30)	BT-PD1-1230(+CS)	12F	.038 inch (0.97mm)
BT-PD1-0635(+CS)	6F	.035 inch (0.89mm)	34 cm (XX=35)	BT-PD1-1235(+CS)	12F	.038 inch (0.97mm)
BT-PD1-0640(+CS)	6F	.035 inch (0.89mm)	39 cm (XX=40)	BT-PD1-1240(+CS)	12F	.038 inch (0.97mm)
BT-PD1-0645(+CS)	6F	.035 inch (0.89mm)	44 cm (XX=45)	BT-PD1-1245(+CS)	12F	.038 inch (0.97mm)
BT-PD1-0650(+CS)	6F	.035 inch (0.89mm)	49 cm (XX=50)	BT-PD1-1250(+CS)	12F	.038 inch (0.97mm)
BT-PD1-0720(+CS)	7F	.035 inch (0.89mm)	19 cm (XX=20)	BT-PD1-1420(+CS)	14F	.038 inch (0.97mm)
BT-PD1-0725(+CS)	7F	.035 inch (0.89mm)	24 cm (XX=25)	BT-PD1-1425(+CS)	14F	.038 inch (0.97mm)
BT-PD1-0730(+CS)	7F	.035 inch (0.89mm)	29 cm (XX=30)	BT-PD1-1430(+CS)	14F	.038 inch (0.97mm)
BT-PD1-0735(+CS)	7F	.035 inch (0.89mm)	34 cm (XX=35)	BT-PD1-1435(+CS)	14F	.038 inch (0.97mm)
BT-PD1-0740(+CS)	7F	.035 inch (0.89mm)	39 cm (XX=40)	BT-PD1-1440(+CS)	14F	.038 inch (0.97mm)
BT-PD1-0745(+CS)	7F	.035 inch (0.89mm)	44 cm (XX=45)	BT-PD1-1445(+CS)	14F	.038 inch (0.97mm)
BT-PD1-0750(+CS)	7F	.035 inch (0.89mm)	49 cm (XX=50)	BT-PD1-1450(+CS)	14F	.038 inch (0.97mm)
BT-PD1-0820(+CS)	8F	.038 inch (0.97mm)	19 cm (XX=20)	BT-PD1-1620(+CS)	16F	.038 inch (0.97mm)
BT-PD1-0825(+CS)	8F	.038 inch (0.97mm)	24 cm (XX=25)	BT-PD1-1625(+CS)	16F	.038 inch (0.97mm)
BT-PD1-0830(+CS)	8F	.038 inch (0.97mm)	29 cm (XX=30)	BT-PD1-1630(+CS)	16F	.038 inch (0.97mm)
BT-PD1-0835(+CS)	8F	.038 inch (0.97mm)	34 cm (XX=35)	BT-PD1-1635(+CS)	16F	.038 inch (0.97mm)
BT-PD1-0840(+CS)	8F	.038 inch (0.97mm)	39 cm (XX=40)	BT-PD1-1640(+CS)	16F	.038 inch (0.97mm)
BT-PD1-0845(+CS)	8F	.038 inch (0.97mm)	44 cm (XX=45)	BT-PD1-1645(+CS)	16F	.038 inch (0.97mm)
BT-PD1-0850(+CS)	8F	.038 inch (0.97mm)	49 cm (XX=50)	BT-PD1-1650(+CS)	16F	.038 inch (0.97mm)
BT-PD1-0920(+CS)	9F	.038 inch (0.97mm)	19 cm (XX=20)			
BT-PD1-0925(+CS)	9F	.038 inch (0.97mm)	24 cm (XX=25)			
BT-PD1-0930(+CS)	9F	.038 inch (0.97mm)	29 cm (XX=30)			
BT-PD1-0935(+CS)	9F	.038 inch (0.97mm)	34 cm (XX=35)			
BT-PD1-0940(+CS)	9F	.038 inch (0.97mm)	39 cm (XX=40)			
BT-PD1-0945(+CS)	9F	.038 inch (0.97mm)	44 cm (XX=45)			
BT-PD1-0950(+CS)	9F	.038 inch (0.97mm)	49 cm (XX=50)			


INTRODUCER SET




INTRODUCER SET

- A kit basically consists of:Co-axial Dilator, Introducer Needle, Guide Wire. See our detailed specifications listed below.
- Introducer set series are legally approved in the Republic of China (Taiwan) to be manufactured, exported and freely sold.

6Fr-4Fr Co-axial Dilator+Introducer Needle+Guide Wire(0.018"; Nitinol)						
Co-axial Dilator	Order Number	Co-axial Dilator			Accessories	
French Size		Working Length	Marker Band	Sheath Accept Guide Wire Diameter	Dilator Accept Guide Wire Diameter	
6F	BT-PDN-618-RB	18cm	with	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet *Radiopaque Coil Tip 0.018 inch(0.46mm) x80 cm Nitinol Guide Wire
6F	BT-PDN-626-RB	26cm	with	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet *Radiopaque Coil Tip 0.018 inch(0.46mm) x80 cm Nitinol Guide Wire
6F	BT-PDN-618	18cm	without	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet *Radiopaque Coil Tip 0.018 inch(0.46mm) x80 cm Nitinol Guide Wire
6F	BT-PDN-626	26cm	without	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet *Radiopaque Coil Tip 0.018 inch(0.46mm) x80 cm Nitinol Guide Wire




6Fr-4Fr Co-axial Dilator+Introducer Needle+Guide Wire(0.018"; Nitinol)+Guide Wire(0.035"; Stainless Steel)						
Co-axial Dilator	Order Number	Co-axial Dilator			Accessories	
French Size		Working Length	Marker Band	Sheath Accept Guide Wire Diameter	Dilator Accept Guide Wire Diameter	
6F	BT-PDN-618-RB(+GW2)	18cm	with	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet *Radiopaque Coil Tip 0.018 inch(0.46mm) x80 cm Nitinol Guide Wire *3mm J tip 0.035 inch(0.89mm)x90cm Guide Wire
6F	BT-PDN-618-RB(+GW2A)	18cm	with	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet *Radiopaque Coil Tip 0.018 inch(0.46mm) x80 cm Nitinol Guide Wire *3mm J tip 0.035 inch(0.89mm)x150cm Guide Wire
6F	BT-PDN-626-RB(+GW2)	26cm	with	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet *Radiopaque Coil Tip 0.018 inch(0.46mm) x80 cm Nitinol Guide Wire *3mm J tip 0.035 inch(0.89mm)x90cm Guide Wire
6F	BT-PDN-626-RB(+GW2A)	26cm	with	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet *Radiopaque Coil Tip 0.018 inch(0.46mm) x80 cm Nitinol Guide Wire *3mm J tip 0.035 inch(0.89mm)x150cm Guide Wire
6F	BT-PDN-618(+GW2)	18cm	without	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet *Radiopaque Coil Tip 0.018 inch(0.46mm) x80 cm Nitinol Guide Wire *3mm J tip 0.035 inch(0.89mm)x90cm Guide Wire
6F	BT-PDN-618(+GW2A)	18cm	without	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet *Radiopaque Coil Tip 0.018 inch(0.46mm) x80 cm Nitinol Guide Wire *3mm J tip 0.035 inch(0.89mm)x150cm Guide Wire
6F	BT-PDN-626(+GW2)	26cm	without	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet *Radiopaque Coil Tip 0.018 inch(0.46mm) x80 cm Nitinol Guide Wire *3mm J tip 0.035 inch(0.89mm)x90cm Guide Wire
6F	BT-PDN-626(+GW2A)	26cm	without	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet *Radiopaque Coil Tip 0.018 inch(0.46mm) x80 cm Nitinol Guide Wire *3mm J tip 0.035 inch(0.89mm)x150cm Guide Wire




INTRODUCER SET

- A kit basically consists of: Co-axial Dilator, Introducer Needle, Guide Wire. See our detailed specifications listed below.
- Introducer set series are legally approved in the Republic of China (Taiwan) to be manufactured, exported and freely sold.

6Fr-4Fr Co-axial Dilator+Introducer Needle						
Co-axial Dilator	Order Number	Co-axial Dilator			Accessories	
French Size		Working Length	Marker Band	Sheath Accept Guide Wire Diameter	Dilator Accept Guide Wire Diameter	
6F	BT-PDN1-618-RB	18cm	with	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet
6F	BT-PDN1-626-RB	26cm	with	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet
6F	BT-PDN1-618	18cm	without	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet
6F	BT-PDN1-626	26cm	without	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet



6Fr-4Fr Co-axial Dilator+Introducer Needle+Guide Wire(0.035"; Stainless Steel)						
Co-axial Dilator	Order Number	Co-axial Dilator			Accessories	
French Size		Working Length	Marker Band	Sheath Accept Guide Wire Diameter	Dilator Accept Guide Wire Diameter	
6F	BT-PDN1-618-RB(+GW2)	18cm	with	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet *3mm J tip 0.035 inch(0.89mm)x90cm Guide Wire
6F	BT-PDN1-618-RB(+GW2A)	18cm	with	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet *3mm J tip 0.035 inch(0.89mm)x150cm Guide Wire
6F	BT-PDN1-626-RB(+GW2)	26cm	with	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet *3mm J tip 0.035 inch(0.89mm)x90cm Guide Wire
6F	BT-PDN1-626-RB(+GW2A)	26cm	with	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet *3mm J tip 0.035 inch(0.89mm)x150cm Guide Wire
6F	BT-PDN1-618(+GW2)	18cm	without	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet *3mm J tip 0.035 inch(0.89mm)x90cm Guide Wire
6F	BT-PDN1-618(+GW2A)	18cm	without	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet *3mm J tip 0.035 inch(0.89mm)x150cm Guide Wire
6F	BT-PDN1-626(+GW2)	26cm	without	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet *3mm J tip 0.035 inch(0.89mm)x90cm Guide Wire
6F	BT-PDN1-626(+GW2A)	26cm	without	0.018" + 0.038"	0.018"	21G(0.82mm) Introducer Needle with stylet *3mm J tip 0.035 inch(0.89mm)x150cm Guide Wire




DRAINAGE CATHETER KIT



INTRODUCER DRAINAGE CATHETER KIT (NK2 Type)

- A kit basically consists of: Pigtail Drainage Catheter with Straightener, Metal Stiffening Cannula, Flexible Stiffening Cannula (only for 7~16F catheter), 6Fr Co-axial Dilator, Dilators, 21G Introducer Needle, 0.018"Guide Wire (80cm), 0.035" Guidewire (90cm) . See our detailed specifications listed below.

NK2 Kit



With SAFETY STRING LOCK						
Order Number	Size	Included Dilators	Working Lengh	Order Number	Size	Included Dilators
BT-PDS-0620-W-NK2-T(+DB/CS)	6F	7F	19 cm (XX=20)	BT-PDS-1020-W-NK2-T(+DB/CS)	10F	8F, 10F, 11F
BT-PDS-0625-W-NK2-T(+DB/CS)	6F	7F	24 cm (XX=25)	BT-PDS-1025-W-NK2-T(+DB/CS)	10F	8F, 10F, 11F
BT-PDS-0630-W-NK2-T(+DB/CS)	6F	7F	29 cm (XX=30)	BT-PDS-1030-W-NK2-T(+DB/CS)	10F	8F, 10F, 11F
BT-PDS-0635-W-NK2-T(+DB/CS)	6F	7F	34 cm (XX=35)	BT-PDS-1035-W-NK2-T(+DB/CS)	10F	8F, 10F, 11F
BT-PDS-0640-W-NK2-T(+DB/CS)	6F	7F	39 cm (XX=40)	BT-PDS-1040-W-NK2-T(+DB/CS)	10F	8F, 10F, 11F
BT-PDS-0645-W-NK2-T(+DB/CS)	6F	7F	44 cm (XX=45)	BT-PDS-1045-W-NK2-T(+DB/CS)	10F	8F, 10F, 11F
BT-PDS-0650-W-NK2-T(+DB/CS)	6F	7F	49 cm (XX=50)	BT-PDS-1050-W-NK2-T(+DB/CS)	10F	8F, 10F, 11F
BT-PDS-0720-W-NK2-T(+DB/CS)	7F	7F, 8F	19 cm (XX=20)	BT-PDS-1220-W-NK2-T(+DB/CS)	12F	8F, 10F, 12F, 13F
BT-PDS-0725-W-NK2-T(+DB/CS)	7F	7F, 8F	24 cm (XX=25)	BT-PDS-1225-W-NK2-T(+DB/CS)	12F	8F, 10F, 12F, 13F
BT-PDS-0730-W-NK2-T(+DB/CS)	7F	7F, 8F	29 cm (XX=30)	BT-PDS-1230-W-NK2-T(+DB/CS)	12F	8F, 10F, 12F, 13F
BT-PDS-0735-W-NK2-T(+DB/CS)	7F	7F, 8F	34 cm (XX=35)	BT-PDS-1235-W-NK2-T(+DB/CS)	12F	8F, 10F, 12F, 13F
BT-PDS-0740-W-NK2-T(+DB/CS)	7F	7F, 8F	39 cm (XX=40)	BT-PDS-1240-W-NK2-T(+DB/CS)	12F	8F, 10F, 12F, 13F
BT-PDS-0745-W-NK2-T(+DB/CS)	7F	7F, 8F	44 cm (XX=45)	BT-PDS-1245-W-NK2-T(+DB/CS)	12F	8F, 10F, 12F, 13F
BT-PDS-0750-W-NK2-T(+DB/CS)	7F	7F, 8F	49 cm (XX=50)	BT-PDS-1250-W-NK2-T(+DB/CS)	12F	8F, 10F, 12F, 13F
BT-PDS-0820-W-NK2-T(+DB/CS)	8F	8F, 9F	19 cm (XX=20)	BT-PDS-1420-W-NK2-T(+DB/CS)	14F	8F, 10F, 12F, 14F, 15F
BT-PDS-0825-W-NK2-T(+DB/CS)	8F	8F, 9F	24 cm (XX=25)	BT-PDS-1425-W-NK2-T(+DB/CS)	14F	8F, 10F, 12F, 14F, 15F
BT-PDS-0830-W-NK2-T(+DB/CS)	8F	8F, 9F	29 cm (XX=30)	BT-PDS-1430-W-NK2-T(+DB/CS)	14F	8F, 10F, 12F, 14F, 15F
BT-PDS-0835-W-NK2-T(+DB/CS)	8F	8F, 9F	34 cm (XX=35)	BT-PDS-1435-W-NK2-T(+DB/CS)	14F	8F, 10F, 12F, 14F, 15F
BT-PDS-0840-W-NK2-T(+DB/CS)	8F	8F, 9F	39 cm (XX=40)	BT-PDS-1440-W-NK2-T(+DB/CS)	14F	8F, 10F, 12F, 14F, 15F
BT-PDS-0845-W-NK2-T(+DB/CS)	8F	8F, 9F	44 cm (XX=45)	BT-PDS-1445-W-NK2-T(+DB/CS)	14F	8F, 10F, 12F, 14F, 15F
BT-PDS-0850-W-NK2-T(+DB/CS)	8F	8F, 9F	49 cm (XX=50)	BT-PDS-1450-W-NK2-T(+DB/CS)	14F	8F, 10F, 12F, 14F, 15F
BT-PDS-0920-W-NK2-T(+DB/CS)	9F	8F, 10F	19 cm (XX=20)	BT-PDS-1620-W-NK2-T(+DB/CS)	16F	8F, 10F, 12F, 14F, 16F
BT-PDS-0925-W-NK2-T(+DB/CS)	9F	8F, 10F	24 cm (XX=25)	BT-PDS-1625-W-NK2-T(+DB/CS)	16F	8F, 10F, 12F, 14F, 16F
BT-PDS-0930-W-NK2-T(+DB/CS)	9F	8F, 10F	29 cm (XX=30)	BT-PDS-1630-W-NK2-T(+DB/CS)	16F	8F, 10F, 12F, 14F, 16F
BT-PDS-0935-W-NK2-T(+DB/CS)	9F	8F, 10F	34 cm (XX=35)	BT-PDS-1635-W-NK2-T(+DB/CS)	16F	8F, 10F, 12F, 14F, 16F
BT-PDS-0940-W-NK2-T(+DB/CS)	9F	8F, 10F	39 cm (XX=40)	BT-PDS-1640-W-NK2-T(+DB/CS)	16F	8F, 10F, 12F, 14F, 16F
BT-PDS-0945-W-NK2-T(+DB/CS)	9F	8F, 10F	44 cm (XX=45)	BT-PDS-1645-W-NK2-T(+DB/CS)	16F	8F, 10F, 12F, 14F, 16F
BT-PDS-0950-W-NK2-T(+DB/CS)	9F	8F, 10F	49 cm (XX=50)	BT-PDS-1650-W-NK2-T(+DB/CS)	16F	8F, 10F, 12F, 14F, 16F

Without SAFETY STRING LOCK

Order Number	Size	Included Dilators	Working Lengh	Order Number	Size	Included Dilators
BT-PDS-0520-NK2-T(+DB/CS)	5F		19 cm (XX=20)	BT-PDS-1020-NK2-T(+DB/CS)	10F	8F, 10F, 11F
BT-PDS-0525-NK2-T(+DB/CS)	5F		24 cm (XX=25)	BT-PDS-1025-NK2-T(+DB/CS)	10F	8F, 10F, 11F
BT-PDS-0530-NK2-T(+DB/CS)	5F		29 cm (XX=30)	BT-PDS-1030-NK2-T(+DB/CS)	10F	8F, 10F, 11F
BT-PDS-0535-NK2-T(+DB/CS)	5F		34 cm (XX=35)	BT-PDS-1035-NK2-T(+DB/CS)	10F	8F, 10F, 11F
BT-PDS-0540-NK2-T(+DB/CS)	5F		39 cm (XX=40)	BT-PDS-1040-NK2-T(+DB/CS)	10F	8F, 10F, 11F
BT-PDS-0545-NK2-T(+DB/CS)	5F		44 cm (XX=45)	BT-PDS-1045-NK2-T(+DB/CS)	10F	8F, 10F, 11F
BT-PDS-0550-NK2-T(+DB/CS)	5F		49 cm (XX=50)	BT-PDS-1050-NK2-T(+DB/CS)	10F	8F, 10F, 11F
BT-PDS-0620-NK2-T(+DB/CS)	6F	7F	19 cm (XX=20)	BT-PDS-1220-NK2-T(+DB/CS)	12F	8F, 10F, 12F, 13F
BT-PDS-0625-NK2-T(+DB/CS)	6F	7F	24 cm (XX=25)	BT-PDS-1225-NK2-T(+DB/CS)	12F	8F, 10F, 12F, 13F
BT-PDS-0630-NK2-T(+DB/CS)	6F	7F	29 cm (XX=30)	BT-PDS-1230-NK2-T(+DB/CS)	12F	8F, 10F, 12F, 13F
BT-PDS-0635-NK2-T(+DB/CS)	6F	7F	34 cm (XX=35)	BT-PDS-1235-NK2-T(+DB/CS)	12F	8F, 10F, 12F, 13F
BT-PDS-0640-NK2-T(+DB/CS)	6F	7F	39 cm (XX=40)	BT-PDS-1240-NK2-T(+DB/CS)	12F	8F, 10F, 12F, 13F
BT-PDS-0645-NK2-T(+DB/CS)	6F	7F	44 cm (XX=45)	BT-PDS-1245-NK2-T(+DB/CS)	12F	8F, 10F, 12F, 13F
BT-PDS-0650-NK2-T(+DB/CS)	6F	7F	49 cm (XX=50)	BT-PDS-1250-NK2-T(+DB/CS)	12F	8F, 10F, 12F, 13F
BT-PDS-0720-NK2-T(+DB/CS)	7F	7F, 8F	19 cm (XX=20)	BT-PDS-1420-NK2-T(+DB/CS)	14F	8F, 10F, 12F, 14F, 15F
BT-PDS-0725-NK2-T(+DB/CS)	7F	7F, 8F	24 cm (XX=25)	BT-PDS-1425-NK2-T(+DB/CS)	14F	8F, 10F, 12F, 14F, 15F
BT-PDS-0730-NK2-T(+DB/CS)	7F	7F, 8F	29 cm (XX=30)	BT-PDS-1430-NK2-T(+DB/CS)	14F	8F, 10F, 12F, 14F, 15F
BT-PDS-0735-NK2-T(+DB/CS)	7F	7F, 8F	34 cm (XX=35)	BT-PDS-1435-NK2-T(+DB/CS)	14F	8F, 10F, 12F, 14F, 15F
BT-PDS-0740-NK2-T(+DB/CS)	7F	7F, 8F	39 cm (XX=40)	BT-PDS-1440-NK2-T(+DB/CS)	14F	8F, 10F, 12F, 14F, 15F
BT-PDS-0745-NK2-T(+DB/CS)	7F	7F, 8F	44 cm (XX=45)	BT-PDS-1445-NK2-T(+DB/CS)	14F	8F, 10F, 12F, 14F, 15F
BT-PDS-0750-NK2-T(+DB/CS)	7F	7F, 8F	49 cm (XX=50)	BT-PDS-1450-NK2-T(+DB/CS)	14F	8F, 10F, 12F, 14F, 15F
BT-PDS-0820-NK2-T(+DB/CS)	8F	8F, 9F	19 cm (XX=20)	BT-PDS-1620-NK2-T(+DB/CS)	16F	8F, 10F, 12F, 14F, 16F
BT-PDS-0825-NK2-T(+DB/CS)	8F	8F, 9F	24 cm (XX=25)	BT-PDS-1625-NK2-T(+DB/CS)	16F	8F, 10F, 12F, 14F, 16F
BT-PDS-0830-NK2-T(+DB/CS)	8F	8F, 9F	29 cm (XX=30)	BT-PDS-1630-NK2-T(+DB/CS)	16F	8F, 10F, 12F, 14F, 16F
BT-PDS-0835-NK2-T(+DB/CS)	8F	8F, 9F	34 cm (XX=35)	BT-PDS-1635-NK2-T(+DB/CS)	16F	8F, 10F, 12F, 14F, 16F
BT-PDS-0840-NK2-T(+DB/CS)	8F	8F, 9F	39 cm (XX=40)	BT-PDS-1640-NK2-T(+DB/CS)	16F	8F, 10F, 12F, 14F, 16F
BT-PDS-0845-NK2-T(+DB/CS)	8F	8F, 9F	44 cm (XX=45)	BT-PDS-1645-NK2-T(+DB/CS)	16F	8F, 10F, 12F, 14F, 16F
BT-PDS-0850-NK2-T(+DB/CS)	8F	8F, 9F	49 cm (XX=50)	BT-PDS-1650-NK2-T(+DB/CS)	16F	8F, 10F, 12F, 14F, 16F
BT-PDS-0920-NK2-T(+DB/CS)	9F	8F, 10F	19 cm (XX=20)			
BT-PDS-0925-NK2-T(+DB/CS)	9F	8F, 10F	24 cm (XX=25)			
BT-PDS-0930-NK2-T(+DB/CS)	9F	8F, 10F	29 cm (XX=30)			
BT-PDS-0935-NK2-T(+DB/CS)	9F	8F, 10F	34 cm (XX=35)			
BT-PDS-0940-NK2-T(+DB/CS)	9F	8F, 10F	39 cm (XX=40)			
BT-PDS-0945-NK2-T(+DB/CS)	9F	8F, 10F	44 cm (XX=45)			
BT-PDS-0950-NK2-T(+DB/CS)	9F	8F, 10F	49 cm (XX=50)			