


UP-6EN

TWS Blue Tooth Hearing Aid

6EN1

6EN2



Various Accessories

Microphone Transmitter



Charging Case



TWS Hearing Aid

Feature

- Fashionable Earphones design.
- Come with a Bluetooth Microphone Transmitter serving as a remote control.
- The transmitter enables clear audio reception for conferences and TV entertainment..
- Handle wider range of frequency and provide high resolution voice and music playing.
- Bluetooth streaming connection.

Specification

- Max Gain: 33 ± 3 dB
- Power Supply: USB Type C
- Continuous Use: 6 hours
- Hearing Aid : 23.5x 19 x 37 mm
- Charge Box: 68 x 29.5 x 37.5 mm
- Weight: 4.5 ± 0.5 g/pc,
Charging Box: 28 ± 5 g/pc
Micro Transmitter: 30 ± 5 g/pc
- Material: Aid, Charger-ABS / Earbud-Silicon
- Color: Black/White

 Recommend for people with
Mild to Moderate hearing loss

Distributed by:
BIOTEQ®