


UP-6S31

ITC Hearing Aid



Invisible and Discreet Design




Fit for either ear



Noise Reduction Design

New and Comfort Experience




Protect Your Ear From
Sudden Burst Sound

Feature

- In-The-Canal
- Long-lasting battery life
- Invisible and discreet design
- Easy to use and comfortable fit
- Easy to remove with a tiny pulling cord

Specification

- Max Gain: 30 ± 5 dB
- Power Supply: 10 Zinc-Air Battery
- Battery Life: 60 hours
- Dimensions: 18.1 x 14 x 11.3 mm
- Weight: 1.5g
- Material: Aid, Case-ABS /
Earbud-Silicone
- Color: Beige

 Recommend for people with
Mild to Moderate hearing loss

Distributed by:
BIOTEQ®