



## **INTRAVENOUS SETS**

## **BIOTEQ®** Intravenous sets meet your needs

We offer a broad range of Intravenous sets to meet your needs in hemodialysis center, general hospitals and other areas.

### **Continuous and precise infusion of liquids and parenteral solutions are important.**

- Precision flow rate control is attainable
- Ease of flow rate maintenance


#### **BI-136 Series for Hemodialysis Centers**

##### **BI-136C**



Infusion set, vented, universal spike:20 drops/ml, roller clamp, male luer lock, soft PVC, approximate set length:140cm, latex free.

##### **BI-136T**




Infusion set, non-vented, universal spike:20 drops/ml, roller clamp, male luer lock, soft PVC, approximate set length:140cm, latex free.





## BI-137 Series For General Hospitals

### BI-137B




20 drops/ml

Infusion set, vented, universal spike:20 drops/ml, roller clamp, male luer connector, soft PVC, DEHP free, approximate set length:190cm, latex free.

### BI-137BY

#### Needleless




20 drops/ml

Infusion set, vented, universal spike:20 drops/ml, roller clamp, needleless Y connector, male luer connector, soft PVC, DEHP free, approximate set length:190cm, latex free.

### BI-137BFY


#### Needleless



20 drops/ml

Infusion set, vented, universal spike:20 drops/ml, roller clamp, needleless Y connector, male luer connector, filter, soft PVC, DEHP free, approximate set length:190cm, latex free.

### BI-137BFD




60 drops/ml

Infusion set, vented, universal spike:60 drops/ml, roller clamp, male luer connector, filter, soft PVC, DEHP free, approximate set length:190cm, latex free.

### BI-137BDY

#### Needleless




60 drops/ml

Infusion set, vented, universal spike:60 drops/ml, roller clamp, needleless Y connector, male luer connector, soft PVC, DEHP free, approximate set length:190cm, latex free.

### BI-137BFDY

#### Needleless




60 drops/ml

Infusion set, vented, universal spike:60 drops/ml, roller clamp, needleless Y connector, male luer connector, filter, soft PVC, DEHP free, approximate set length:190cm, latex free.

## BI-137 Series For General Hospitals

### BI-137




20 drops/ml

Infusion set, vented, universal spike:20 drops/ml, roller clamp, male luer connector, soft PVC, approximate set length:190cm, latex free.

### BI-137Y

#### Needleless




20 drops/ml

Infusion set, vented, universal spike:20 drops/ml, roller clamp, needleless Y connector, male luer connector, soft PVC, approximate set length:190cm, latex free.

### BI-137FY


#### Needleless



20 drops/ml

Infusion set, vented, universal spike:20 drops/ml, roller clamp, needleless Y connector, male luer connector, filter, soft PVC, approximate set length:190cm, latex free.

### BI-137FD




60 drops/ml

Infusion set, vented, universal spike:60 drops/ml, roller clamp, male luer lock connector, filter, soft PVC, approximate set length:190cm, latex free.

### BI-137DY

#### Needleless




60 drops/ml

Infusion set, vented, universal spike:60 drops/ml, roller clamp, needleless Y connector, male luer connector, soft PVC, approximate set length:190cm, latex free.

### BI-137FDY

#### Needleless




60 drops/ml

Infusion set, vented, universal spike:60 drops/ml, roller clamp, needleless Y connector, male luer connector, filter, soft PVC, approximate set length:190cm, latex free.

## BI-137 Series For General Hospitals


### BI-137NY



20 drops/ml

Infusion set, vented universal spike: 20 drops/ml, roller clamp, Y connector, male luer connector, soft PVC, approximate set length: 190cm, latex free


### BI-137NFY



20 drops/ml

Infusion set, vented universal spike: 20 drops/ml, roller clamp, Y connector, male luer connector, filter, soft PVC, approximate set length: 190cm, latex free


### BI-137NDY



60 drops/ml

Infusion set, vented universal spike: 60 drops/ml, roller clamp, Y connector, male luer connector, soft PVC, approximate set length: 190cm, latex free

### BI-137NFDY




60 drops/ml

Infusion set, vented universal spike: 60 drops/ml, roller clamp, Y connector, male luer connector, filter, soft PVC, approximate set length: 190cm, latex free

## BT-163 Series : 100ml Precision Infusion Set


### BT-163FN1

Needleless



Precision Infusion set, 100 ml precision chamber (assembled needleless connector upper lid), microdrip chamber 60 drops/ml, spike 20 drops/ml, air vented spike, roller clamp, male luer connector, needleless Y connector, soft PVC, DEHP free, approximate set length: 200cm, latex free.

### BT-163FY1



Precision Infusion set, 100 ml precision chamber, microdrip chamber 60 drops/ml, spike 20 drops/ml, air vented spike, roller clamp, male luer connector, rubber Y connector, soft PVC, DEHP free, approximate set length: 200cm, latex free