




IV Solution Bag

- Standard soft container for pharmaceutical IV solution
- Custom-made printing
- Superior high resistance to impact during transportation
- No air intake, free of particle contamination, micro-organisms and pyrogens
- Quality assurance in compliance with ISO 9001, ISO 13485, CE0434 and GMP
- Economy of cost-effectiveness



Twist off




Stopper

Capacity in ml:

50/100/150/200/250/500/1000/1500/2000/3000/5000

BIOTEQ IV Infusion bag

- Range of Capacity: 50cc to 5,000cc.
- Range of Length: 125mm to 400mm.
- Interval (Pitch) between stopper and twist-off is 25mm & 40mm.
- Tube length of stopper (F) and twist-off (E) can be custom made according customer's requirements.



IV Infusion Bag Specifications:

Item	Capacity	Width(A)	Length (B)	Thickness(C)	Pitch (D)
1	50	80	125	0.32	40
2	100	80	143	0.32	40
3	150	90	180	0.35	40
4	150	102	152.5	0.32	40
5	200	102	162	0.32	40
6	250	120	165	0.32	40
7	500	120	192	0.32	40
8	500	120	210	0.32	40
9	500	120	220	0.32	40
10	500	120	225	0.32	40
11	1,000	120	314	0.32	40
12	1,000	120	325	0.32	40
13	1,000	127	307.5	0.35	40
14	1,500	127	400	0.35	40
15	2,000	180	335	0.35	25
16	2,000	180	345	0.35	25
17	3,000	210	340	0.35	25
18	3,000	210	380	0.35	25
19	3,000	210	390	0.35	25
20	5,000	270	380	0.35	25



For more information, please contact:
sales@bioteq.com.tw

BIOTEQUE CORPORATION
5F-6, No.23, Chang-An E. Rd.,
Taipei 104, Taiwan
Phone: (886)-2-2571-0269, Fax: (886)-2-2536-1967
Http://www.bioteq.com.tw

ISO 13485, GMP, CE certified
@Bioteque Corp. April, 2010
Brochure number : BT-IVB-01