

# BIOTEQ<sup>®</sup>

## HSG Catheter



**For Hysterosalpingography  
or Sonohysterography  
procedures**



### Features and Benefits

1. Two lumen configuration: One lumen is for contrast saline, the other is for balloon inflation.
2. Made of 100% silicone material, soft, smooth and comfortable, non-toxic medical silica gel material that do not produce a toxic or immune response within the human body..
3. The balloon wall is thin, firm and smooth ensuring the proper seal of cervical os.
4. Introducer tube can be shaped to facilitate fitting the contour of the body.
3. Open tip, Closed tip with sideport are available.

