



DRAINAGE CATHETER SET-DIRECT ACCESS



- Used for nephrostomy, abscess, and other multipurpose drainage.
- Catheter placement is performed using one step (direct access) technique.
- The safety string lock system provides strong fixation of the pigtail.
- Hydrophilic decreased friction of the distal catheter facilitates easy insertion.
- Large oval drainage holes increase drainage capability.
- Supplied sterile in peel-open packages.
- Intended for one-time use.

PIGTAIL TYPE (BT-PD1-series) + FLEXIBLE STIFFENING CANNULA(FSC)



With SAFETY STRING LOCK						
Order Number	Size	Accepts Guide Wire Dia.	Working Lengh	Order Number	Size	Accepts Guide Wire Dia.
BT-PD1-0720-W(+FSC)	7F	.035 inch (0.89mm)	19 cm (XX=20)	BT-PD1-1220-W(+FSC)	12F	.038 inch (0.97mm)
BT-PD1-0725-W(+FSC)	7F	.035 inch (0.89mm)	24 cm (XX=25)	BT-PD1-1225-W(+FSC)	12F	.038 inch (0.97mm)
BT-PD1-0730-W(+FSC)	7F	.035 inch (0.89mm)	29 cm (XX=30)	BT-PD1-1230-W(+FSC)	12F	.038 inch (0.97mm)
BT-PD1-0735-W(+FSC)	7F	.035 inch (0.89mm)	34 cm (XX=35)	BT-PD1-1235-W(+FSC)	12F	.038 inch (0.97mm)
BT-PD1-0740-W(+FSC)	7F	.035 inch (0.89mm)	39 cm (XX=40)	BT-PD1-1240-W(+FSC)	12F	.038 inch (0.97mm)
BT-PD1-0745-W(+FSC)	7F	.035 inch (0.89mm)	44 cm (XX=45)	BT-PD1-1245-W(+FSC)	12F	.038 inch (0.97mm)
BT-PD1-0750-W(+FSC)	7F	.035 inch (0.89mm)	49 cm (XX=50)	BT-PD1-1250-W(+FSC)	12F	.038 inch (0.97mm)
BT-PD1-0820-W(+FSC)	8F	.038 inch (0.97mm)	19 cm (XX=20)	BT-PD1-1420-W(+FSC)	14F	.038 inch (0.97mm)
BT-PD1-0825-W(+FSC)	8F	.038 inch (0.97mm)	24 cm (XX=25)	BT-PD1-1425-W(+FSC)	14F	.038 inch (0.97mm)
BT-PD1-0830-W(+FSC)	8F	.038 inch (0.97mm)	29 cm (XX=30)	BT-PD1-1430-W(+FSC)	14F	.038 inch (0.97mm)
BT-PD1-0835-W(+FSC)	8F	.038 inch (0.97mm)	34 cm (XX=35)	BT-PD1-1435-W(+FSC)	14F	.038 inch (0.97mm)
BT-PD1-0840-W(+FSC)	8F	.038 inch (0.97mm)	39 cm (XX=40)	BT-PD1-1440-W(+FSC)	14F	.038 inch (0.97mm)
BT-PD1-0845-W(+FSC)	8F	.038 inch (0.97mm)	44 cm (XX=45)	BT-PD1-1445-W(+FSC)	14F	.038 inch (0.97mm)
BT-PD1-0850-W(+FSC)	8F	.038 inch (0.97mm)	49 cm (XX=50)	BT-PD1-1450-W(+FSC)	14F	.038 inch (0.97mm)
BT-PD1-0920-W(+FSC)	9F	.038 inch (0.97mm)	19 cm (XX=20)	BT-PD1-1620-W(+FSC)	16F	.038 inch (0.97mm)
BT-PD1-0925-W(+FSC)	9F	.038 inch (0.97mm)	24 cm (XX=25)	BT-PD1-1625-W(+FSC)	16F	.038 inch (0.97mm)
BT-PD1-0930-W(+FSC)	9F	.038 inch (0.97mm)	29 cm (XX=30)	BT-PD1-1630-W(+FSC)	16F	.038 inch (0.97mm)
BT-PD1-0935-W(+FSC)	9F	.038 inch (0.97mm)	34 cm (XX=35)	BT-PD1-1635-W(+FSC)	16F	.038 inch (0.97mm)
BT-PD1-0940-W(+FSC)	9F	.038 inch (0.97mm)	39 cm (XX=40)	BT-PD1-1640-W(+FSC)	16F	.038 inch (0.97mm)
BT-PD1-0945-W(+FSC)	9F	.038 inch (0.97mm)	44 cm (XX=45)	BT-PD1-1645-W(+FSC)	16F	.038 inch (0.97mm)
BT-PD1-0950-W(+FSC)	9F	.038 inch (0.97mm)	49 cm (XX=50)	BT-PD1-1650-W(+FSC)	16F	.038 inch (0.97mm)
BT-PD1-1020-W(+FSC)	10F	.038 inch (0.97mm)	19 cm (XX=20)			
BT-PD1-1025-W(+FSC)	10F	.038 inch (0.97mm)	24 cm (XX=25)			
BT-PD1-1030-W(+FSC)	10F	.038 inch (0.97mm)	29 cm (XX=30)			
BT-PD1-1035-W(+FSC)	10F	.038 inch (0.97mm)	34 cm (XX=35)			
BT-PD1-1040-W(+FSC)	10F	.038 inch (0.97mm)	39 cm (XX=40)			
BT-PD1-1045-W(+FSC)	10F	.038 inch (0.97mm)	44 cm (XX=45)			
BT-PD1-1050-W(+FSC)	10F	.038 inch (0.97mm)	49 cm (XX=50)			

Without SAFETY STRING LOCK						
Order Number	Size	Accepts Guide Wire Dia.	Working Lengh	Order Number	Size	Accepts Guide Wire Dia.
BT-PD1-0720(+FSC)	7F	.035 inch (0.89mm)	19 cm (XX=20)	BT-PD1-1220(+FSC)	12F	.038 inch (0.97mm)
BT-PD1-0725(+FSC)	7F	.035 inch (0.89mm)	24 cm (XX=25)	BT-PD1-1225(+FSC)	12F	.038 inch (0.97mm)
BT-PD1-0730(+FSC)	7F	.035 inch (0.89mm)	29 cm (XX=30)	BT-PD1-1230(+FSC)	12F	.038 inch (0.97mm)
BT-PD1-0735(+FSC)	7F	.035 inch (0.89mm)	34 cm (XX=35)	BT-PD1-1235(+FSC)	12F	.038 inch (0.97mm)
BT-PD1-0740(+FSC)	7F	.035 inch (0.89mm)	39 cm (XX=40)	BT-PD1-1240(+FSC)	12F	.038 inch (0.97mm)
BT-PD1-0745(+FSC)	7F	.035 inch (0.89mm)	44 cm (XX=45)	BT-PD1-1245(+FSC)	12F	.038 inch (0.97mm)
BT-PD1-0750(+FSC)	7F	.035 inch (0.89mm)	49 cm (XX=50)	BT-PD1-1250(+FSC)	12F	.038 inch (0.97mm)
BT-PD1-0820(+FSC)	8F	.038 inch (0.97mm)	19 cm (XX=20)	BT-PD1-1420(+FSC)	14F	.038 inch (0.97mm)
BT-PD1-0825(+FSC)	8F	.038 inch (0.97mm)	24 cm (XX=25)	BT-PD1-1425(+FSC)	14F	.038 inch (0.97mm)
BT-PD1-0830(+FSC)	8F	.038 inch (0.97mm)	29 cm (XX=30)	BT-PD1-1430(+FSC)	14F	.038 inch (0.97mm)
BT-PD1-0835(+FSC)	8F	.038 inch (0.97mm)	34 cm (XX=35)	BT-PD1-1435(+FSC)	14F	.038 inch (0.97mm)
BT-PD1-0840(+FSC)	8F	.038 inch (0.97mm)	39 cm (XX=40)	BT-PD1-1440(+FSC)	14F	.038 inch (0.97mm)
BT-PD1-0845(+FSC)	8F	.038 inch (0.97mm)	44 cm (XX=45)	BT-PD1-1445(+FSC)	14F	.038 inch (0.97mm)
BT-PD1-0850(+FSC)	8F	.038 inch (0.97mm)	49 cm (XX=50)	BT-PD1-1450(+FSC)	14F	.038 inch (0.97mm)
BT-PD1-0920(+FSC)	9F	.038 inch (0.97mm)	19 cm (XX=20)	BT-PD1-1620(+FSC)	14F	.038 inch (0.97mm)
BT-PD1-0925(+FSC)	9F	.038 inch (0.97mm)	24 cm (XX=25)	BT-PD1-1625(+FSC)	16F	.038 inch (0.97mm)
BT-PD1-0930(+FSC)	9F	.038 inch (0.97mm)	29 cm (XX=30)	BT-PD1-1630(+FSC)	16F	.038 inch (0.97mm)
BT-PD1-0935(+FSC)	9F	.038 inch (0.97mm)	34 cm (XX=35)	BT-PD1-1635(+FSC)	16F	.038 inch (0.97mm)
BT-PD1-0940(+FSC)	9F	.038 inch (0.97mm)	39 cm (XX=40)	BT-PD1-1640(+FSC)	16F	.038 inch (0.97mm)
BT-PD1-0945(+FSC)	9F	.038 inch (0.97mm)	44 cm (XX=45)	BT-PD1-1645(+FSC)	16F	.038 inch (0.97mm)
BT-PD1-0950(+FSC)	9F	.038 inch (0.97mm)	49 cm (XX=50)	BT-PD1-1650(+FSC)	16F	.038 inch (0.97mm)
BT-PD1-1020(+FSC)	10F	.038 inch (0.97mm)	19 cm (XX=20)			
BT-PD1-1025(+FSC)	10F	.038 inch (0.97mm)	24 cm (XX=25)			
BT-PD1-1030(+FSC)	10F	.038 inch (0.97mm)	29 cm (XX=30)			
BT-PD1-1035(+FSC)	10F	.038 inch (0.97mm)	34 cm (XX=35)			
BT-PD1-1040(+FSC)	10F	.038 inch (0.97mm)	39 cm (XX=40)			
BT-PD1-1045(+FSC)	10F	.038 inch (0.97mm)	44 cm (XX=45)			
BT-PD1-1050(+FSC)	10F	.038 inch (0.97mm)	49 cm (XX=50)			