



DRAINAGE CATHETER SET-DIRECT ACCESS



- Used for abscess, and other multipurpose drainage.
- Catheter placement is performed using one step (direct access) technique.
- Hydrophilic coating decreased friction of the distal catheter facilitates easy insertion.
- Large oval drainage holes increase drainage capability.
- Supplied sterile in peel-open packages.
- Intended for one-time use.

CURVED (J) TYPE (BT-PD1-series(J))



Without SAFETY STRING LOCK

Order Number	Size	Accepts Guide Wire Dia.	Working Lengh	Order Number	Size	Accepts Guide Wire Dia.
BT-PD1-0520(J)	5F	.035 inch (0.89mm)	19 cm (XX=20)	BT-PD1-1020(J)	10F	.038 inch (0.97mm)
BT-PD1-0525(J)	5F	.035 inch (0.89mm)	24 cm (XX=25)	BT-PD1-1025(J)	10F	.038 inch (0.97mm)
BT-PD1-0530(J)	5F	.035 inch (0.89mm)	29 cm (XX=30)	BT-PD1-1030(J)	10F	.038 inch (0.97mm)
BT-PD1-0535(J)	5F	.035 inch (0.89mm)	34 cm (XX=35)	BT-PD1-1035(J)	10F	.038 inch (0.97mm)
BT-PD1-0540(J)	5F	.035 inch (0.89mm)	39 cm (XX=40)	BT-PD1-1040(J)	10F	.038 inch (0.97mm)
BT-PD1-0545(J)	5F	.035 inch (0.89mm)	44 cm (XX=45)	BT-PD1-1045(J)	10F	.038 inch (0.97mm)
BT-PD1-0550(J)	5F	.035 inch (0.89mm)	49 cm (XX=50)	BT-PD1-1050(J)	10F	.038 inch (0.97mm)
BT-PD1-0620(J)	6F	.035 inch (0.89mm)	19 cm (XX=20)	BT-PD1-1220(J)	12F	.038 inch (0.97mm)
BT-PD1-0625(J)	6F	.035 inch (0.89mm)	24 cm (XX=25)	BT-PD1-1225(J)	12F	.038 inch (0.97mm)
BT-PD1-0630(J)	6F	.035 inch (0.89mm)	29 cm (XX=30)	BT-PD1-1230(J)	12F	.038 inch (0.97mm)
BT-PD1-0635(J)	6F	.035 inch (0.89mm)	34 cm (XX=35)	BT-PD1-1235(J)	12F	.038 inch (0.97mm)
BT-PD1-0640(J)	6F	.035 inch (0.89mm)	39 cm (XX=40)	BT-PD1-1240(J)	12F	.038 inch (0.97mm)
BT-PD1-0645(J)	6F	.035 inch (0.89mm)	44 cm (XX=45)	BT-PD1-1245(J)	12F	.038 inch (0.97mm)
BT-PD1-0650(J)	6F	.035 inch (0.89mm)	49 cm (XX=50)	BT-PD1-1250(J)	12F	.038 inch (0.97mm)
BT-PD1-0720(J)	7F	.035 inch (0.89mm)	19 cm (XX=20)	BT-PD1-1420(J)	14F	.038 inch (0.97mm)
BT-PD1-0725(J)	7F	.035 inch (0.89mm)	24 cm (XX=25)	BT-PD1-1425(J)	14F	.038 inch (0.97mm)
BT-PD1-0730(J)	7F	.035 inch (0.89mm)	29 cm (XX=30)	BT-PD1-1430(J)	14F	.038 inch (0.97mm)
BT-PD1-0735(J)	7F	.035 inch (0.89mm)	34 cm (XX=35)	BT-PD1-1435(J)	14F	.038 inch (0.97mm)
BT-PD1-0740(J)	7F	.035 inch (0.89mm)	39 cm (XX=40)	BT-PD1-1440(J)	14F	.038 inch (0.97mm)
BT-PD1-0745(J)	7F	.035 inch (0.89mm)	44 cm (XX=45)	BT-PD1-1445(J)	14F	.038 inch (0.97mm)
BT-PD1-0750(J)	7F	.035 inch (0.89mm)	49 cm (XX=50)	BT-PD1-1450(J)	14F	.038 inch (0.97mm)
BT-PD1-0820(J)	8F	.038 inch (0.97mm)	19 cm (XX=20)	BT-PD1-1620(J)	16F	.038 inch (0.97mm)
BT-PD1-0825(J)	8F	.038 inch (0.97mm)	24 cm (XX=25)	BT-PD1-1625(J)	16F	.038 inch (0.97mm)
BT-PD1-0830(J)	8F	.038 inch (0.97mm)	29 cm (XX=30)	BT-PD1-1630(J)	16F	.038 inch (0.97mm)
BT-PD1-0835(J)	8F	.038 inch (0.97mm)	34 cm (XX=35)	BT-PD1-1635(J)	16F	.038 inch (0.97mm)
BT-PD1-0840(J)	8F	.038 inch (0.97mm)	39 cm (XX=40)	BT-PD1-1640(J)	16F	.038 inch (0.97mm)
BT-PD1-0845(J)	8F	.038 inch (0.97mm)	44 cm (XX=45)	BT-PD1-1645(J)	16F	.038 inch (0.97mm)
BT-PD1-0850(J)	8F	.038 inch (0.97mm)	49 cm (XX=50)	BT-PD1-1650(J)	16F	.038 inch (0.97mm)
BT-PD1-0920(J)	9F	.038 inch (0.97mm)	19 cm (XX=20)			
BT-PD1-0925(J)	9F	.038 inch (0.97mm)	24 cm (XX=25)			
BT-PD1-0930(J)	9F	.038 inch (0.97mm)	29 cm (XX=30)			
BT-PD1-0935(J)	9F	.038 inch (0.97mm)	34 cm (XX=35)			
BT-PD1-0940(J)	9F	.038 inch (0.97mm)	39 cm (XX=40)			
BT-PD1-0945(J)	9F	.038 inch (0.97mm)	44 cm (XX=45)			
BT-PD1-0950(J)	9F	.038 inch (0.97mm)	49 cm (XX=50)			