


DRAINAGE CATHETER SET-SELDINGER



- Used for abscess, and other multipurpose drainage.
- Catheter placement is performed using Seldinger technique.
- Hydrophilic coating decreased friction of the distal catheter facilitates easy insertion.
- Large oval drainage holes increase drainage capability.
- Supplied sterile in peel-open packages.
- Intended for one-time use.

STRAIGHT TYPE (BT-PD2-series)



Without SAFETY STRING LOCK

Order Number	Size	Accepts Guide Wire Dia.	Working Length	Order Number	Size	Accepts Guide Wire Dia.
BT-PD2-0520	5F	.035 inch (0.89mm)	19 cm (XX=20)	BT-PD2-1020	10F	.038 inch (0.97mm)
BT-PD2-0525	5F	.035 inch (0.89mm)	24 cm (XX=25)	BT-PD2-1025	10F	.038 inch (0.97mm)
BT-PD2-0530	5F	.035 inch (0.89mm)	29 cm (XX=30)	BT-PD2-1030	10F	.038 inch (0.97mm)
BT-PD2-0535	5F	.035 inch (0.89mm)	34 cm (XX=35)	BT-PD2-1035	10F	.038 inch (0.97mm)
BT-PD2-0540	5F	.035 inch (0.89mm)	39 cm (XX=40)	BT-PD2-1040	10F	.038 inch (0.97mm)
BT-PD2-0545	5F	.035 inch (0.89mm)	44 cm (XX=45)	BT-PD2-1045	10F	.038 inch (0.97mm)
BT-PD2-0550	5F	.035 inch (0.89mm)	49 cm (XX=50)	BT-PD2-1050	10F	.038 inch (0.97mm)
BT-PD2-0620	6F	.035 inch (0.89mm)	19 cm (XX=20)	BT-PD2-1220	12F	.038 inch (0.97mm)
BT-PD2-0625	6F	.035 inch (0.89mm)	24 cm (XX=25)	BT-PD2-1225	12F	.038 inch (0.97mm)
BT-PD2-0630	6F	.035 inch (0.89mm)	29 cm (XX=30)	BT-PD2-1230	12F	.038 inch (0.97mm)
BT-PD2-0635	6F	.035 inch (0.89mm)	34 cm (XX=35)	BT-PD2-1235	12F	.038 inch (0.97mm)
BT-PD2-0640	6F	.035 inch (0.89mm)	39 cm (XX=40)	BT-PD2-1240	12F	.038 inch (0.97mm)
BT-PD2-0645	6F	.035 inch (0.89mm)	44 cm (XX=45)	BT-PD2-1245	12F	.038 inch (0.97mm)
BT-PD2-0650	6F	.035 inch (0.89mm)	49 cm (XX=50)	BT-PD2-1250	12F	.038 inch (0.97mm)
BT-PD2-0720	7F	.035 inch (0.89mm)	19 cm (XX=20)	BT-PD2-1420	14F	.038 inch (0.97mm)
BT-PD2-0725	7F	.035 inch (0.89mm)	24 cm (XX=25)	BT-PD2-1425	14F	.038 inch (0.97mm)
BT-PD2-0730	7F	.035 inch (0.89mm)	29 cm (XX=30)	BT-PD2-1430	14F	.038 inch (0.97mm)
BT-PD2-0735	7F	.035 inch (0.89mm)	34 cm (XX=35)	BT-PD2-1435	14F	.038 inch (0.97mm)
BT-PD2-0740	7F	.035 inch (0.89mm)	39 cm (XX=40)	BT-PD2-1440	14F	.038 inch (0.97mm)
BT-PD2-0745	7F	.035 inch (0.89mm)	44 cm (XX=45)	BT-PD2-1445	14F	.038 inch (0.97mm)
BT-PD2-0750	7F	.035 inch (0.89mm)	49 cm (XX=50)	BT-PD2-1450	14F	.038 inch (0.97mm)
BT-PD2-0820	8F	.038 inch (0.97mm)	19 cm (XX=20)	BT-PD2-1620	16F	.038 inch (0.97mm)
BT-PD2-0825	8F	.038 inch (0.97mm)	24 cm (XX=25)	BT-PD2-1625	16F	.038 inch (0.97mm)
BT-PD2-0830	8F	.038 inch (0.97mm)	29 cm (XX=30)	BT-PD2-1630	16F	.038 inch (0.97mm)
BT-PD2-0835	8F	.038 inch (0.97mm)	34 cm (XX=35)	BT-PD2-1635	16F	.038 inch (0.97mm)
BT-PD2-0840	8F	.038 inch (0.97mm)	39 cm (XX=40)	BT-PD2-1640	16F	.038 inch (0.97mm)
BT-PD2-0845	8F	.038 inch (0.97mm)	44 cm (XX=45)	BT-PD2-1645	16F	.038 inch (0.97mm)
BT-PD2-0850	8F	.038 inch (0.97mm)	49 cm (XX=50)	BT-PD2-1650	16F	.038 inch (0.97mm)
BT-PD2-0920	9F	.038 inch (0.97mm)	19 cm (XX=20)			
BT-PD2-0925	9F	.038 inch (0.97mm)	24 cm (XX=25)			
BT-PD2-0930	9F	.038 inch (0.97mm)	29 cm (XX=30)			
BT-PD2-0935	9F	.038 inch (0.97mm)	34 cm (XX=35)			
BT-PD2-0940	9F	.038 inch (0.97mm)	39 cm (XX=40)			
BT-PD2-0945	9F	.038 inch (0.97mm)	44 cm (XX=45)			
BT-PD2-0950	9F	.038 inch (0.97mm)	49 cm (XX=50)			