



# DRAINAGE CATHETER SET-SELDINGER



- Used for abscess, and other multipurpose drainage.
- Catheter placement is performed using Seldinger technique.
- Hydrophilic coating decreased friction of the distal catheter facilitates easy insertion.
- Supplied sterile in peel-open packages.
- Intended for one-time use.

## MALECOT TYPE (BT-PDS-series(M))



### Without SAFETY STRING LOCK

Order Number	Size	Accepts Guide Wire Dia.	Working Length	Order Number	Size	Accepts Guide Wire Dia.
BT-PDS-0620(M)	6F	.035 inch (0.89mm)	19 cm (XX=20)	BT-PDS-1020(M)	10F	.038 inch (0.97mm)
BT-PDS-0625(M)	6F	.035 inch (0.89mm)	24 cm (XX=25)	BT-PDS-1025(M)	10F	.038 inch (0.97mm)
BT-PDS-0630(M)	6F	.035 inch (0.89mm)	29 cm (XX=30)	BT-PDS-1030(M)	10F	.038 inch (0.97mm)
BT-PDS-0635(M)	6F	.035 inch (0.89mm)	34 cm (XX=35)	BT-PDS-1035(M)	10F	.038 inch (0.97mm)
BT-PDS-0640(M)	6F	.035 inch (0.89mm)	39 cm (XX=40)	BT-PDS-1040(M)	10F	.038 inch (0.97mm)
BT-PDS-0645(M)	6F	.035 inch (0.89mm)	44 cm (XX=45)	BT-PDS-1045(M)	10F	.038 inch (0.97mm)
BT-PDS-0650(M)	6F	.035 inch (0.89mm)	49 cm (XX=50)	BT-PDS-1050(M)	10F	.038 inch (0.97mm)
BT-PDS-0720(M)	7F	.035 inch (0.89mm)	19 cm (XX=20)	BT-PDS-1220(M)	12F	.038 inch (0.97mm)
BT-PDS-0725(M)	7F	.035 inch (0.89mm)	24 cm (XX=25)	BT-PDS-1225(M)	12F	.038 inch (0.97mm)
BT-PDS-0730(M)	7F	.035 inch (0.89mm)	29 cm (XX=30)	BT-PDS-1230(M)	12F	.038 inch (0.97mm)
BT-PDS-0735(M)	7F	.035 inch (0.89mm)	34 cm (XX=35)	BT-PDS-1235(M)	12F	.038 inch (0.97mm)
BT-PDS-0740(M)	7F	.035 inch (0.89mm)	39 cm (XX=40)	BT-PDS-1240(M)	12F	.038 inch (0.97mm)
BT-PDS-0745(M)	7F	.035 inch (0.89mm)	44 cm (XX=45)	BT-PDS-1245(M)	12F	.038 inch (0.97mm)
BT-PDS-0750(M)	7F	.035 inch (0.89mm)	49 cm (XX=50)	BT-PDS-1250(M)	12F	.038 inch (0.97mm)
BT-PDS-0820(M)	8F	.038 inch (0.97mm)	19 cm (XX=20)	BT-PDS-1420(M)	14F	.038 inch (0.97mm)
BT-PDS-0825(M)	8F	.038 inch (0.97mm)	24 cm (XX=25)	BT-PDS-1425(M)	14F	.038 inch (0.97mm)
BT-PDS-0830(M)	8F	.038 inch (0.97mm)	29 cm (XX=30)	BT-PDS-1430(M)	14F	.038 inch (0.97mm)
BT-PDS-0835(M)	8F	.038 inch (0.97mm)	34 cm (XX=35)	BT-PDS-1435(M)	14F	.038 inch (0.97mm)
BT-PDS-0840(M)	8F	.038 inch (0.97mm)	39 cm (XX=40)	BT-PDS-1440(M)	14F	.038 inch (0.97mm)
BT-PDS-0845(M)	8F	.038 inch (0.97mm)	44 cm (XX=45)	BT-PDS-1445(M)	14F	.038 inch (0.97mm)
BT-PDS-0850(M)	8F	.038 inch (0.97mm)	49 cm (XX=50)	BT-PDS-1450(M)	14F	.038 inch (0.97mm)
BT-PDS-0920(M)	9F	.038 inch (0.97mm)	19 cm (XX=20)	BT-PDS-1620(M)	16F	.038 inch (0.97mm)
BT-PDS-0925(M)	9F	.038 inch (0.97mm)	24 cm (XX=25)	BT-PDS-1625(M)	16F	.038 inch (0.97mm)
BT-PDS-0930(M)	9F	.038 inch (0.97mm)	29 cm (XX=30)	BT-PDS-1630(M)	16F	.038 inch (0.97mm)
BT-PDS-0935(M)	9F	.038 inch (0.97mm)	34 cm (XX=35)	BT-PDS-1635(M)	16F	.038 inch (0.97mm)
BT-PDS-0940(M)	9F	.038 inch (0.97mm)	39 cm (XX=40)	BT-PDS-1640(M)	16F	.038 inch (0.97mm)
BT-PDS-0945(M)	9F	.038 inch (0.97mm)	44 cm (XX=45)	BT-PDS-1645(M)	16F	.038 inch (0.97mm)
BT-PDS-0950(M)	9F	.038 inch (0.97mm)	49 cm (XX=50)	BT-PDS-1650(M)	16F	.038 inch (0.97mm)