


DRAINAGE CATHETER SET-SELDINGER



- Used for nephrostomy, abscess, and other multipurpose drainage.
- Catheter placement is performed using Seldinger technique.
- The safety string lock system provides strong fixation of the pigtail.
- Hydrophilic coating decreased friction of the distal catheter facilitates easy insertion.
- Large oval drainage holes increase drainage capability.
- Radiopaque band aids in proper placement of the catheter.
- Supplied sterile in peel-open packages.

CLOSE LOOP PIGTAIL TYPE (BT-PDS-series(CL)-RB)



With SAFETY STRING LOCK						
Order Number	Size	Accepts Guide Wire Dia.	Working Length	Order Number	Size	Accepts Guide Wire Dia.
BT-PDS-0620-W(CL)-RB	6F	.035 inch (0.89mm)	19 cm (XX=20)	BT-PDS-1220-W(CL)-RB	12F	.038 inch (0.97mm)
BT-PDS-0625-W(CL)-RB	6F	.035 inch (0.89mm)	24 cm (XX=25)	BT-PDS-1225-W(CL)-RB	12F	.038 inch (0.97mm)
BT-PDS-0630-W(CL)-RB	6F	.035 inch (0.89mm)	29 cm (XX=30)	BT-PDS-1230-W(CL)-RB	12F	.038 inch (0.97mm)
BT-PDS-0635-W(CL)-RB	6F	.035 inch (0.89mm)	34 cm (XX=35)	BT-PDS-1235-W(CL)-RB	12F	.038 inch (0.97mm)
BT-PDS-0640-W(CL)-RB	6F	.035 inch (0.89mm)	39 cm (XX=40)	BT-PDS-1240-W(CL)-RB	12F	.038 inch (0.97mm)
BT-PDS-0645-W(CL)-RB	6F	.035 inch (0.89mm)	44 cm (XX=45)	BT-PDS-1245-W(CL)-RB	12F	.038 inch (0.97mm)
BT-PDS-0650-W(CL)-RB	6F	.035 inch (0.89mm)	49 cm (XX=50)	BT-PDS-1250-W(CL)-RB	12F	.038 inch (0.97mm)
BT-PDS-0820-W(CL)-RB	8F	.038 inch (0.97mm)	19 cm (XX=20)	BT-PDS-1420-W(CL)-RB	14F	.038 inch (0.97mm)
BT-PDS-0825-W(CL)-RB	8F	.038 inch (0.97mm)	24 cm (XX=25)	BT-PDS-1425-W(CL)-RB	14F	.038 inch (0.97mm)
BT-PDS-0830-W(CL)-RB	8F	.038 inch (0.97mm)	29 cm (XX=30)	BT-PDS-1430-W(CL)-RB	14F	.038 inch (0.97mm)
BT-PDS-0835-W(CL)-RB	8F	.038 inch (0.97mm)	34 cm (XX=35)	BT-PDS-1435-W(CL)-RB	14F	.038 inch (0.97mm)
BT-PDS-0840-W(CL)-RB	8F	.038 inch (0.97mm)	39 cm (XX=40)	BT-PDS-1440-W(CL)-RB	14F	.038 inch (0.97mm)
BT-PDS-0845-W(CL)-RB	8F	.038 inch (0.97mm)	44 cm (XX=45)	BT-PDS-1445-W(CL)-RB	14F	.038 inch (0.97mm)
BT-PDS-0850-W(CL)-RB	8F	.038 inch (0.97mm)	49 cm (XX=50)	BT-PDS-1450-W(CL)-RB	14F	.038 inch (0.97mm)
BT-PDS-1020-W(CL)-RB	10F	.038 inch (0.97mm)	19 cm (XX=20)			
BT-PDS-1025-W(CL)-RB	10F	.038 inch (0.97mm)	24 cm (XX=25)			
BT-PDS-1030-W(CL)-RB	10F	.038 inch (0.97mm)	29 cm (XX=30)			
BT-PDS-1035-W(CL)-RB	10F	.038 inch (0.97mm)	34 cm (XX=35)			
BT-PDS-1040-W(CL)-RB	10F	.038 inch (0.97mm)	39 cm (XX=40)			
BT-PDS-1045-W(CL)-RB	10F	.038 inch (0.97mm)	44 cm (XX=45)			
BT-PDS-1050-W(CL)-RB	10F	.038 inch (0.97mm)	49 cm (XX=50)			

Without SAFETY STRING LOCK						
Order Number	Size	Accepts Guide Wire Dia.	Working Length	Order Number	Size	Accepts Guide Wire Dia.
BT-PDS-0620(CL)-RB	6F	.035 inch (0.89mm)	19 cm (XX=20)	BT-PDS-1220(CL)-RB	12F	.038 inch (0.97mm)
BT-PDS-0625(CL)-RB	6F	.035 inch (0.89mm)	24 cm (XX=25)	BT-PDS-1225(CL)-RB	12F	.038 inch (0.97mm)
BT-PDS-0630(CL)-RB	6F	.035 inch (0.89mm)	29 cm (XX=30)	BT-PDS-1230(CL)-RB	12F	.038 inch (0.97mm)
BT-PDS-0635(CL)-RB	6F	.035 inch (0.89mm)	34 cm (XX=35)	BT-PDS-1235(CL)-RB	12F	.038 inch (0.97mm)
BT-PDS-0640(CL)-RB	6F	.035 inch (0.89mm)	39 cm (XX=40)	BT-PDS-1240(CL)-RB	12F	.038 inch (0.97mm)
BT-PDS-0645(CL)-RB	6F	.035 inch (0.89mm)	44 cm (XX=45)	BT-PDS-1245(CL)-RB	12F	.038 inch (0.97mm)
BT-PDS-0650(CL)-RB	6F	.035 inch (0.89mm)	49 cm (XX=50)	BT-PDS-1250(CL)-RB	12F	.038 inch (0.97mm)
BT-PDS-0820(CL)-RB	8F	.038 inch (0.97mm)	19 cm (XX=20)	BT-PDS-1420(CL)-RB	14F	.038 inch (0.97mm)
BT-PDS-0825(CL)-RB	8F	.038 inch (0.97mm)	24 cm (XX=25)	BT-PDS-1425(CL)-RB	14F	.038 inch (0.97mm)
BT-PDS-0830(CL)-RB	8F	.038 inch (0.97mm)	29 cm (XX=30)	BT-PDS-1430(CL)-RB	14F	.038 inch (0.97mm)
BT-PDS-0835(CL)-RB	8F	.038 inch (0.97mm)	34 cm (XX=35)	BT-PDS-1435(CL)-RB	14F	.038 inch (0.97mm)
BT-PDS-0840(CL)-RB	8F	.038 inch (0.97mm)	39 cm (XX=40)	BT-PDS-1440(CL)-RB	14F	.038 inch (0.97mm)
BT-PDS-0845(CL)-RB	8F	.038 inch (0.97mm)	44 cm (XX=45)	BT-PDS-1445(CL)-RB	14F	.038 inch (0.97mm)
BT-PDS-0850(CL)-RB	8F	.038 inch (0.97mm)	49 cm (XX=50)	BT-PDS-1450(CL)-RB	14F	.038 inch (0.97mm)
BT-PDS-1020(CL)-RB	10F	.038 inch (0.97mm)	19 cm (XX=20)			
BT-PDS-1025(CL)-RB	10F	.038 inch (0.97mm)	24 cm (XX=25)			
BT-PDS-1030(CL)-RB	10F	.038 inch (0.97mm)	29 cm (XX=30)			
BT-PDS-1035(CL)-RB	10F	.038 inch (0.97mm)	34 cm (XX=35)			
BT-PDS-1040(CL)-RB	10F	.038 inch (0.97mm)	39 cm (XX=40)			
BT-PDS-1045(CL)-RB	10F	.038 inch (0.97mm)	44 cm (XX=45)			
BT-PDS-1050(CL)-RB	10F	.038 inch (0.97mm)	49 cm (XX=50)			