


PERCUTANEOUS DRAINAGE- ACCESSORIES

- Supplied sterile in peel-open packages.
- Intended for one-time use.

DRAINAGE BAG




Order Number	Product Name	Capacity	Description
LU-500-DF	LEG BAG	500ml	Male Luer Lock
LU-800-DF	PCN BAG	800ml	Male Luer Lock
LU-800(F)-DF	PCN BAG	800ml	Female Luer Lock
LU-800(W/FTF-90)-DF	PCN BAG	800ml	Male Luer Lock;With Female/Female Luer Lock Adapter
DB-800-DF	DRAINAGE BAG	800ml	Female Luer Lock Tubing Length Adjustable
DB-800C-MS-V-DF	DRAINAGE BAG	800ml	Male Luer Lock;With V.:Elastic Cloth Belt (Velcvo) S.:Stopcock C.:Back cloth
UB-121T-DF	DISPOSABLE DRAINAGE BAG	2000ml	T-Valve; Olive
UB-121T(FL)-DF	DISPOSABLE DRAINAGE BAG	2000ml	T-Valve; Female Luer Lock
UB-121T(ML)-DF	DISPOSABLE DRAINAGE BAG	2000ml	T-Valve; Male Luer Lock

PERCUTANEOUS DRAINAGE- ACCESSORIES

GUIDEWIRE

- Intended for one-time use.



Nitinol Guidewire



Teflon Coated Guidewire

Order Number	Description
BT-GW-S-018-80	0.018"/80cm/Straight Tip/Palladium Radiopaque Coil Tip/Nitinol Shaft
BT-GW-S-018-80S	0.018"/80cm/Straight Tip/Palladium Radiopaque Coil Tip/Stainless Steel Shaft
BT-GW-J3-035-90	0.035"/90cm/3mm J/Stiff shaft/Teflon Coated
BT-GW-J3-035-90A	0.035"/90cm/3mm J/Teflon Coated
BT-GW-J3-035-150	0.035"/150cm/3mm J/Stiff shaft/Teflon Coated
BT-GW-J3-035-150A	0.035"/150cm/3mm J/Teflon Coated
BT-GW-J3-038-90	0.038"/90cm/3mm J/Teflon Coated
BT-GW-J3-038-150	0.038"/150cm/3mm J/Teflon Coated
BT-GW-S-035-90-HC	0.035"/90cm/Straight/Nitinol/Hydrophilic Coated
BT-GW-S-035-150-HC	0.035"/150cm/Straight/Nitinol/Hydrophilic Coated
BT-GW-S-038-90-HC	0.038"/90cm/Straight/Nitinol/Hydrophilic Coated
BT-GW-S-038-150-HC	0.038"/150cm/Straight/Nitinol/Hydrophilic Coated

SCALPEL

- Supplied sterile in peel-open packages.
- Intended for one-time use.



Order Number	Description
BT-SCL-BN11	No. 11/Minor