



# PERCUTANEOUS DRAINAGE- ACCESSORIES

## 2 PART DILATOR

- Intended for one-time use.
- Dilator material:
- ER:Polyethylene with Radiopaque.
- UPR:Outer layer/Polyurethane with Radiopaque, Inner layer/Polypropylene with Radiopaque.



2-PART DILATORS				
French Size	Order Number	Needle/Dilator Accept Guide Wire Diameter	Sheath Accept Guide Wire Diameter	Working Length
6F-4F	DR-UPR-2P-0610	Dilators:0.018 inch(0.46mm)	0.038 inch(0.97mm)	10
5F-3F	DR-UPR-2P-0510	Dilators:0.018 inch(0.46mm)	0.038 inch(0.97mm)	10
5F-3F	DR-ER-2P-0518	Needle:0.018 inch(0.46mm)	0.038 inch(0.97mm)	18
5F-3F	DR-ER-2P-0526	Needle:0.018 inch(0.46mm)	0.038 inch(0.97mm)	26
4.5F-3F	DR-UPR-2P-4510	Dilators:0.018 inch(0.46mm)	0.038 inch(0.97mm)	10
4F-3F	DR-UPR-2P-0410	Dilators:0.018 inch(0.46mm)	0.038 inch(0.97mm)	10