


PERCUTANEOUS DRAINAGE- ACCESSORIES

3 PART DILATOR

- Intended for one-time use.
- Dilator material:ne with Radiopaque.
- UER:Outer layer/Polyurethane with Radiopaque, Inner layer/Polyethylene with Radiopaque.



3-PART DILATORS

French Size	Order Number	Needle/Dilator Accept Guide Wire Diameter	Sheath Accept Guide Wire Diameter	Working Length
6F-4F	DR-UER-3P-618	Needle:0.018 inch(0.46mm)	0.18"(0.46mm)+0.038"(0.97mm)	18
6F-4F	DR-UER-3P-626	Needle:0.018 inch(0.46mm)	0.18"(0.46mm)+0.038"(0.97mm)	26
6F-4F	DR-UER-3P-0618-RB	Needle:0.018 inch(0.46mm)	0.18"(0.46mm)+0.038"(0.97mm)	18
6F-4F	DR-UER-3P-0626-RB	Needle:0.018 inch(0.46mm)	0.18"(0.46mm)+0.038"(0.97mm)	26