


INTRODUCER DRAINAGE CATHETER KIT (NK1 Type)

- A kit basically consists of: Biliary Catheter with Radiopaque Band with Straightener, Metal Stiffening Cannula, Flexible Stiffening Cannula (only for 7~16F catheter), Dilators, 18G Introducer Needle, 0.035" Guide Wire (150cm). See our detailed specifications listed below.

NK1 Kit



With SAFETY STRING LOCK						
Order Number	Size	Included Dilators	Working Lengh	Order Number	Size	Included Dilators
BT-PDS-0620-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	19 cm (XX=20)	BT-PDS-1220-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0625-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	24 cm (XX=25)	BT-PDS-1225-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0630-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	29 cm (XX=30)	BT-PDS-1230-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0635-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	34 cm (XX=35)	BT-PDS-1235-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0640-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	39 cm (XX=40)	BT-PDS-1240-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0645-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	6F	6F,7F	44 cm (XX=45)	BT-PDS-1245-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0650-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	49 cm (XX=50)	BT-PDS-1250-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 12F, 13F
BT-PDS-0820-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	19 cm (XX=20)	BT-PDS-1420-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0825-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	24 cm (XX=25)	BT-PDS-1425-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0830-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	29 cm (XX=30)	BT-PDS-1430-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0835-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	34 cm (XX=35)	BT-PDS-1435-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0840-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	39 cm (XX=40)	BT-PDS-1440-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0845-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	44 cm (XX=45)	BT-PDS-1445-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0850-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	49 cm (XX=50)	BT-PDS-1450-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-1020-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F	19 cm (XX=20)			
BT-PDS-1025-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F	24 cm (XX=25)			
BT-PDS-1030-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F	29 cm (XX=30)			
BT-PDS-1035-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F	34 cm (XX=35)			
BT-PDS-1040-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F	39 cm (XX=40)			
BT-PDS-1045-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F	44 cm (XX=45)			
BT-PDS-1050-W(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F	49 cm (XX=50)			

Without SAFETY STRING LOCK

Order Number	Size	Included Dilators	Working Lengh	Order Number	Size	Included Dilators
BT-PDS-0620-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	19 cm (XX=20)	BT-PDS-1220-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0625-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	24 cm (XX=25)	BT-PDS-1225-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0630-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	29 cm (XX=30)	BT-PDS-1230-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0635-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	34 cm (XX=35)	BT-PDS-1235-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0640-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	39 cm (XX=40)	BT-PDS-1240-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0645-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	6F	6F,7F	44 cm (XX=45)	BT-PDS-1245-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0650-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	49 cm (XX=50)	BT-PDS-1250-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 12F, 13F
BT-PDS-0820-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	19 cm (XX=20)	BT-PDS-1420-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0825-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	24 cm (XX=25)	BT-PDS-1425-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0830-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	29 cm (XX=30)	BT-PDS-1430-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0835-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	34 cm (XX=35)	BT-PDS-1435-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0840-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	39 cm (XX=40)	BT-PDS-1440-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0845-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	44 cm (XX=45)	BT-PDS-1445-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0850-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	49 cm (XX=50)	BT-PDS-1450-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-1020-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F	19 cm (XX=20)			
BT-PDS-1025-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F	24 cm (XX=25)			
BT-PDS-1030-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F	29 cm (XX=30)			
BT-PDS-1035-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F	34 cm (XX=35)			
BT-PDS-1040-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F	39 cm (XX=40)			
BT-PDS-1045-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F	44 cm (XX=45)			
BT-PDS-1050-(B)-RB-NK1(L)-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F	49 cm (XX=50)			