


INTRODUCER DRAINAGE CATHETER KIT (NK2 Type)

- A kit basically consists of: Pigtail Drainage Catheter with Straightener, Metal Stiffening Cannula, Flexible Stiffening Cannula (only for 7~16F catheter), 6Fr Co-axial Dilator, Dilators, 21G Introducer Needle, 0.018"Guide Wire (80cm), 0.035" Guidewire (90cm) . See our detailed specifications listed below.

NK2 Kit



With SAFETY STRING LOCK						
Order Number	Size	Included Dilators	Working Lengh	Order Number	Size	Included Dilators
BT-PDS-0620-W-NK2-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	19 cm (XX=20)	BT-PDS-1020-W-NK2-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F
BT-PDS-0625-W-NK2-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	24 cm (XX=25)	BT-PDS-1025-W-NK2-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F
BT-PDS-0630-W-NK2-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	29 cm (XX=30)	BT-PDS-1030-W-NK2-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F
BT-PDS-0635-W-NK2-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	34 cm (XX=35)	BT-PDS-1035-W-NK2-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F
BT-PDS-0640-W-NK2-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	39 cm (XX=40)	BT-PDS-1040-W-NK2-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F
BT-PDS-0645-W-NK2-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	44 cm (XX=45)	BT-PDS-1045-W-NK2-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F
BT-PDS-0650-W-NK2-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	49 cm (XX=50)	BT-PDS-1050-W-NK2-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F
BT-PDS-0720-W-NK2-T(+CT/DB/SC/ST/CS/CN)	7F	6F, 7F, 8F	19 cm (XX=20)	BT-PDS-1220-W-NK2-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0725-W-NK2-T(+CT/DB/SC/ST/CS/CN)	7F	6F, 7F, 8F	24 cm (XX=25)	BT-PDS-1225-W-NK2-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0730-W-NK2-T(+CT/DB/SC/ST/CS/CN)	7F	6F, 7F, 8F	29 cm (XX=30)	BT-PDS-1230-W-NK2-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0735-W-NK2-T(+CT/DB/SC/ST/CS/CN)	7F	6F, 7F, 8F	34 cm (XX=35)	BT-PDS-1235-W-NK2-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0740-W-NK2-T(+CT/DB/SC/ST/CS/CN)	7F	6F, 7F, 8F	39 cm (XX=40)	BT-PDS-1240-W-NK2-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0745-W-NK2-T(+CT/DB/SC/ST/CS/CN)	7F	6F, 7F, 8F	44 cm (XX=45)	BT-PDS-1245-W-NK2-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0750-W-NK2-T(+CT/DB/SC/ST/CS/CN)	7F	6F, 8F, 9F	49 cm (XX=50)	BT-PDS-1250-W-NK2-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0820-W-NK2-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	19 cm (XX=20)	BT-PDS-1420-W-NK2-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0825-W-NK2-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	24 cm (XX=25)	BT-PDS-1425-W-NK2-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0830-W-NK2-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	29 cm (XX=30)	BT-PDS-1430-W-NK2-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0835-W-NK2-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	34 cm (XX=35)	BT-PDS-1435-W-NK2-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0840-W-NK2-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	39 cm (XX=40)	BT-PDS-1440-W-NK2-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0845-W-NK2-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	44 cm (XX=45)	BT-PDS-1445-W-NK2-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0850-W-NK2-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	49 cm (XX=50)	BT-PDS-1450-W-NK2-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0920-W-NK2-T(+CT/DB/SC/ST/CS/CN)	9F	6F, 8F, 10F	19 cm (XX=20)	BT-PDS-1620-W-NK2-T(+CT/DB/SC/ST/CS/CN)	16F	6F, 8F, 10F, 12F, 14F, 16F
BT-PDS-0925-W-NK2-T(+CT/DB/SC/ST/CS/CN)	9F	6F, 8F, 10F	24 cm (XX=25)	BT-PDS-1625-W-NK2-T(+CT/DB/SC/ST/CS/CN)	16F	6F, 8F, 10F, 12F, 14F, 16F
BT-PDS-0930-W-NK2-T(+CT/DB/SC/ST/CS/CN)	9F	6F, 8F, 10F	29 cm (XX=30)	BT-PDS-1630-W-NK2-T(+CT/DB/SC/ST/CS/CN)	16F	6F, 8F, 10F, 12F, 14F, 16F
BT-PDS-0935-W-NK2-T(+CT/DB/SC/ST/CS/CN)	9F	6F, 8F, 10F	34 cm (XX=35)	BT-PDS-1635-W-NK2-T(+CT/DB/SC/ST/CS/CN)	16F	6F, 8F, 10F, 12F, 14F, 16F
BT-PDS-0940-W-NK2-T(+CT/DB/SC/ST/CS/CN)	9F	6F, 8F, 10F	39 cm (XX=40)	BT-PDS-1640-W-NK2-T(+CT/DB/SC/ST/CS/CN)	16F	6F, 8F, 10F, 12F, 14F, 16F
BT-PDS-0945-W-NK2-T(+CT/DB/SC/ST/CS/CN)	9F	6F, 8F, 10F	44 cm (XX=45)	BT-PDS-1645-W-NK2-T(+CT/DB/SC/ST/CS/CN)	16F	6F, 8F, 10F, 12F, 14F, 16F
BT-PDS-0950-W-NK2-T(+CT/DB/SC/ST/CS/CN)	9F	6F, 8F, 10F	49 cm (XX=50)	BT-PDS-1650-W-NK2-T(+CT/DB/SC/ST/CS/CN)	16F	6F, 8F, 10F, 12F, 14F, 16F

Without SAFETY STRING LOCK						
Order Number	Size	Included Dilators	Working Lengh	Order Number	Size	Included Dilators
BT-PDS-0520-NK2-T(+CT/DB/SC/ST/CS/CN)	5F	6F	19 cm (XX=20)	BT-PDS-1020-NK2-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F
BT-PDS-0525-NK2-T(+CT/DB/SC/ST/CS/CN)	5F	6F	24 cm (XX=25)	BT-PDS-1025-NK2-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F
BT-PDS-0530-NK2-T(+CT/DB/SC/ST/CS/CN)	5F	6F	29 cm (XX=30)	BT-PDS-1030-NK2-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F
BT-PDS-0535-NK2-T(+CT/DB/SC/ST/CS/CN)	5F	6F	34 cm (XX=35)	BT-PDS-1035-NK2-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F
BT-PDS-0540-NK2-T(+CT/DB/SC/ST/CS/CN)	5F	6F	39 cm (XX=40)	BT-PDS-1040-NK2-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F
BT-PDS-0545-NK2-T(+CT/DB/SC/ST/CS/CN)	5F	6F	44 cm (XX=45)	BT-PDS-1045-NK2-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F
BT-PDS-0550-NK2-T(+CT/DB/SC/ST/CS/CN)	5F	6F	49 cm (XX=50)	BT-PDS-1050-NK2-T(+CT/DB/SC/ST/CS/CN)	10F	6F, 8F, 10F, 11F
BT-PDS-0620-NK2-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	19 cm (XX=20)	BT-PDS-1220-NK2-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0625-NK2-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	24 cm (XX=25)	BT-PDS-1225-NK2-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0630-NK2-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	29 cm (XX=30)	BT-PDS-1230-NK2-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0635-NK2-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	34 cm (XX=35)	BT-PDS-1235-NK2-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0640-NK2-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	39 cm (XX=40)	BT-PDS-1240-NK2-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0645-NK2-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	44 cm (XX=45)	BT-PDS-1245-NK2-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0650-NK2-T(+CT/DB/SC/ST/CS/CN)	6F	6F, 7F	49 cm (XX=50)	BT-PDS-1250-NK2-T(+CT/DB/SC/ST/CS/CN)	12F	6F, 8F, 10F, 12F, 13F
BT-PDS-0720-NK2-T(+CT/DB/SC/ST/CS/CN)	7F	6F, 7F, 8F	19 cm (XX=20)	BT-PDS-1420-NK2-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0725-NK2-T(+CT/DB/SC/ST/CS/CN)	7F	6F, 7F, 8F	24 cm (XX=25)	BT-PDS-1425-NK2-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0730-NK2-T(+CT/DB/SC/ST/CS/CN)	7F	6F, 7F, 8F	29 cm (XX=30)	BT-PDS-1430-NK2-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0735-NK2-T(+CT/DB/SC/ST/CS/CN)	7F	6F, 7F, 8F	34 cm (XX=35)	BT-PDS-1435-NK2-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0740-NK2-T(+CT/DB/SC/ST/CS/CN)	7F	6F, 7F, 8F	39 cm (XX=40)	BT-PDS-1440-NK2-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0745-NK2-T(+CT/DB/SC/ST/CS/CN)	7F	6F, 7F, 8F	44 cm (XX=45)	BT-PDS-1445-NK2-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0750-NK2-T(+CT/DB/SC/ST/CS/CN)	7F	6F, 7F, 8F	49 cm (XX=50)	BT-PDS-1450-NK2-T(+CT/DB/SC/ST/CS/CN)	14F	6F, 8F, 10F, 12F, 14F, 15F
BT-PDS-0820-NK2-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	19 cm (XX=20)	BT-PDS-1620-NK2-T(+CT/DB/SC/ST/CS/CN)	16F	6F, 8F, 10F, 12F, 14F, 16F
BT-PDS-0825-NK2-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	24 cm (XX=25)	BT-PDS-1625-NK2-T(+CT/DB/SC/ST/CS/CN)	16F	6F, 8F, 10F, 12F, 14F, 16F
BT-PDS-0830-NK2-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	29 cm (XX=30)	BT-PDS-1630-NK2-T(+CT/DB/SC/ST/CS/CN)	16F	6F, 8F, 10F, 12F, 14F, 16F
BT-PDS-0835-NK2-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	34 cm (XX=35)	BT-PDS-1635-NK2-T(+CT/DB/SC/ST/CS/CN)	16F	6F, 8F, 10F, 12F, 14F, 16F
BT-PDS-0840-NK2-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	39 cm (XX=40)	BT-PDS-1640-NK2-T(+CT/DB/SC/ST/CS/CN)	16F	6F, 8F, 10F, 12F, 14F, 16F
BT-PDS-0845-NK2-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	44 cm (XX=45)	BT-PDS-1645-NK2-T(+CT/DB/SC/ST/CS/CN)	16F	6F, 8F, 10F, 12F, 14F, 16F
BT-PDS-0850-NK2-T(+CT/DB/SC/ST/CS/CN)	8F	6F, 8F, 9F	49 cm (XX=50)	BT-PDS-1650-NK2-T(+CT/DB/SC/ST/CS/CN)	16F	6F, 8F, 10F, 12F, 14F, 16F
BT-PDS-0920-NK2-T(+CT/DB/SC/ST/CS/CN)	9F	6F, 8F, 10F	19 cm (XX=20)			
BT-PDS-0925-NK2-T(+CT/DB/SC/ST/CS/CN)	9F	6F, 8F, 10F	24 cm (XX=25)			
BT-PDS-0930-NK2-T(+CT/DB/SC/ST/CS/CN)	9F	6F, 8F, 10F	29 cm (XX=30)			
BT-PDS-0935-NK2-T(+CT/DB/SC/ST/CS/CN)	9F	6F, 8F, 10F	34 cm (XX=35)			
BT-PDS-0940-NK2-T(+CT/DB/SC/ST/CS/CN)	9F	6F, 8F, 10F	39 cm (XX=40)			
BT-PDS-0945-NK2-T(+CT/DB/SC/ST/CS/CN)	9F	6F, 8F, 10F	44 cm (XX=45)			
BT-PDS-0950-NK2-T(+CT/DB/SC/ST/CS/CN)	9F	6F, 8F, 10F	49 cm (XX=50)			