


# 24-Hour Adult T-piece

**BIOTEC®** The 24-Hr adult T-piece has a 90° angle for gas flow exit, it is packaged with a cap. T-piece port can be attached with flex tube, converting the T-piece design into a small-volume reservoir.

A patient can maintain effective oxygenation delivered by the ventilator while being suctioned whether in the hospital or at home.



**Clear T-Piece**

Attachment to the artificial airway with 22 mm connections on both sides.



Attaching Flex tube to one side of T-piece reduces volume of incoming air increase SaO2, connect ventilator breathing circuit to the other side.

**One-way Irrigation Port**


Helps not only maintain the pressure of the ventilator circuit but also takes the injection of normal saline to clean up the catheter after each suction procedure.

**Soft Sleeve**

Highly compact and soft sleeve to effectively prevent infiltration of blood and sputum.



Model	Catheter Size	Length(mm)	Control Button	Type	Packaging
BT-SC1-1054NT	10 Fr/3.3mm	540	Black	T-Piece Endotracheal Open Irrigation Port	10sets/case
BT-SC1-1254NT	12 Fr/4.0mm		White		
BT-SC1-1454NT	14 Fr/4.6mm		Green		
BT-SC1-1654NT	16 Fr/5.3mm		Orange		
BT-SC2-1054NT	10 Fr/3.3mm	540	Black	T-Piece Endotracheal One-way Needleless Irrigation Port	10sets/case
BT-SC2-1254NT	12 Fr/4.0mm		White		
BT-SC2-1454NT	14 Fr/4.6mm		Green		
BT-SC2-1654NT	16 Fr/5.3mm		Orange		
BT-SC7-1054NT	10 Fr/3.3mm	540	Black	T-Piece Endotracheal Luer Lock Irrigation Port	10sets/case
BT-SC7-1254NT	12 Fr/4.0mm		White		
BT-SC7-1454NT	14 Fr/4.6mm		Green		
BT-SC7-1654NT	16 Fr/5.3mm		Orange		
BT-SC1-1036NT	10 Fr/3.3mm	360	Black	T-Piece Tracheostomy Open Irrigation Port	10sets/case
BT-SC1-1236NT	12 Fr/4.0mm		White		
BT-SC1-1436NT	14 Fr/4.6mm		Green		
BT-SC1-1636NT	16 Fr/5.3mm		Orange		
BT-SC2-1036NT	10 Fr/3.3mm	360	Black	T-Piece Tracheostomy One-way Needleless Irrigation Port	10sets/case
BT-SC2-1236NT	12 Fr/4.0mm		White		
BT-SC2-1436NT	14 Fr/4.6mm		Green		
BT-SC2-1636NT	16 Fr/5.3mm		Orange		
BT-SC7-1036NT	10 Fr/3.3mm	360	Black	T-Piece Tracheostomy Luer Lock Irrigation Port	10sets/case
BT-SC7-1236NT	12 Fr/4.0mm		White		
BT-SC7-1436NT	14 Fr/4.6mm		Green		
BT-SC7-1636NT	16 Fr/5.3mm		Orange		



ACCESSORY CODE:(If Flex tube or Disconnecting Wedge is desired, then include a W, F, WF in Model.#)  
BT-SCX-XXXXNTWF:24Hr T-Piece type+W(Disconnecting Wedge)+F(Flex Tube)

REMARK: Products on this page are not CE-marked.