

# 24-Hour Pediatric Y Connector

**BIOTEQ®** 24-Hour Pediatric Catheter set offers pediatric/neonatal intubated patients a more comfortable suction procedure without breaking up the oxygen support. With color marked tubing and softened catheter tip, the depth of catheter insertion can be precisely controlled and the damage on the mucus membrane can be minimized to the lowest level.

**Day-of-change Sticker**  
Reminds caregivers the correct day for changing a new catheter.



**Soft Sleeve**  
Fully covers the suction catheter for better protection on both healthcare practitioners and the patients. Meanwhile, the transparency and the softness of the material offer easy operation for the procedure.

**Y Adapters**  
Ranging from 2.0/2.5/3.0/3.5/4.0/4.5/5.0/5.5mm are specially prepared for pediatric patients.



**Ergonomic Control Button**  
Takes only single move to manipulate the suction source.



**Color-Marked Tubing**  
Provides better recognition for insertion depth.



**Black Number Coded Tubing**  
Provides better recognition for insertion depth.



Model	Catheter Size	Length(mm)	Control Button	Y Adapters	Type	Packaging
BT-SC1-0532NC	05Fr/1.6mm	320	Gray	2.0mm/2.5mm	Color Open Irrigation Port	10sets/case
BT-SC1-0632NC	06Fr/2.0mm		Light Green	2.5mm/3.0mm/3.5mm		
BT-SC1-0732NC	07Fr/2.3mm		Purple	3.0mm/3.5mm/4.0mm		
BT-SC1-0836NC	08Fr/2.6mm		Light Blue	3.5mm/4.0mm/4.5mm		
BT-SC1-1041NC	10Fr/3.3mm	415	Black	4.5mm/5.0mm/5.5mm	Color One-way Needleless Irrigation Port	10sets/case
BT-SC2-0532NC	05Fr/1.6mm	320	Gray	2.0mm/2.5mm		
BT-SC2-0632NC	06Fr/2.0mm		Light Green	2.5mm/3.0mm/3.5mm		
BT-SC2-0732NC	07Fr/2.3mm		Purple	3.0mm/3.5mm/4.0mm		
BT-SC2-0836NC	08Fr/2.6mm	360	Light Blue	3.5mm/4.0mm/4.5mm		
BT-SC2-1041NC	10Fr/3.3mm	415	Black	4.5mm/5.0mm/5.5mm	Color Luer Lock Irrigation Port	10sets/case
BT-SC7-0532NC	05Fr/1.6mm	320	Gray	2.0mm/2.5mm		
BT-SC7-0632NC	06Fr/2.0mm		Light Green	2.5mm/3.0mm/3.5mm		
BT-SC7-0732NC	07Fr/2.3mm		Purple	3.0mm/3.5mm/4.0mm		
BT-SC7-0836NC	08Fr/2.6mm	360	Light Blue	3.5mm/4.0mm/4.5mm		
BT-SC7-1041NC	10Fr/3.3mm	415	Black	4.5mm/5.0mm/5.5mm		



ACCESSORY CODE:(If Flex tube is desired, then include a C, F, CF in Model#.)  
C: Pediatric Color catheter,  
If black number coded graduations catheters are desired, then goes without "C".  
BT-SCX-XXXXNCF:24 Hr Pediatric type+C(Pediatric Color catheter)+F(Flex Tube)

