

Suction Jet 72-Hour Double Swivel Elbow Stop Valve Type

Smooth Tip
Minimizes the damage on mucous membrane to the lowest level.

Double Swivel Connectors
Allow 360 degrees adjustment on tubing angles to improve patients' comfort.



MDI (Metered Dose Inhaler) Port
Allows delivery of aerosol drugs as bronchodilators like albuterol into the ventilator circuit to mechanically ventilated patients



Day-of-change Sticker
Reminds caregivers the correct day for changing a new catheter.



Ergonomic Control Button
Takes only single move to manipulate the suction source.



Special Stop Valve
Stands over 500 times insertion and withdrawal manipulation.



Keep closed chamber for flush cleaning the catheter tip.

When the suction is on, the closed chamber makes strong flush flowing through and perfectly clean the catheter tip.

Luer lock Irrigation Port
Provides needle-free access via a syringe to help prevent needle sticks.



Model	Catheter Size	Length(mm)	Control Button	Type	Packaging
BT-SC1-1054ND7	10 Fr/3.3mm	540	Black	Double Swivel Endotracheal Open Irrigation Port	10sets/case
BT-SC1-1254ND7	12 Fr/4.0mm		White		
BT-SC1-1454ND7	14 Fr/4.6mm		Green		
BT-SC1-1654ND7	16 Fr/5.3mm		Orange		
BT-SC1-1854ND7	18 Fr/6.0mm		Red		
BT-SC2-1054ND7	10 Fr/3.3mm	540	Black	Double Swivel Endotracheal One-way Needleless Irrigation Port	10sets/case
BT-SC2-1254ND7	12 Fr/4.0mm		White		
BT-SC2-1454ND7	14 Fr/4.6mm		Green		
BT-SC2-1654ND7	16 Fr/5.3mm		Orange		
BT-SC2-1854ND7	18 Fr/6.0mm		Red		
BT-SC5-1054ND7	10 Fr/3.3mm	540	Black	Double Swivel Endotracheal One-way Needleless Irrigation Port MDI Port	10sets/case
BT-SC5-1254ND7	12 Fr/4.0mm		White		
BT-SC5-1454ND7	14 Fr/4.6mm		Green		
BT-SC5-1654ND7	16 Fr/5.3mm		Orange		
BT-SC5-1854ND7	18 Fr/6.0mm		Red		
BT-SC7-1054ND7	10 Fr/3.3mm	540	Black	Double Swivel Endotracheal Luer Lock Irrigation Port	10sets/case
BT-SC7-1254ND7	12 Fr/4.0mm		White		
BT-SC7-1454ND7	14 Fr/4.6mm		Green		
BT-SC7-1654ND7	16 Fr/5.3mm		Orange		
BT-SC7-1854ND7	18 Fr/6.0mm		Red		
BT-SC8-1054ND7	10 Fr/3.3mm	540	Black	Double Swivel Endotracheal Luer Lock Irrigation Port MDI Port	10sets/case
BT-SC8-1254ND7	12 Fr/4.0mm		White		
BT-SC8-1454ND7	14 Fr/4.6mm		Green		
BT-SC8-1654ND7	16 Fr/5.3mm		Orange		
BT-SC8-1854ND7	18 Fr/6.0mm		Red		
BT-SC1-1036ND7	10 Fr/3.3mm	360	Black	Double Swivel Tracheostomy Open Irrigation Port	10sets/case
BT-SC1-1236ND7	12 Fr/4.0mm		White		
BT-SC1-1436ND7	14 Fr/4.6mm		Green		
BT-SC1-1636ND7	16 FR/5.3mm		Orange		
BT-SC1-1836ND7	18 Fr/6.0mm		Red		
BT-SC2-1036ND7	10 Fr/3.3mm	360	Black	Double Swivel Tracheostomy One-way Needleless Irrigation Port	10sets/case
BT-SC2-1236ND7	12 Fr/4.0mm		White		
BT-SC2-1436ND7	14 Fr/4.6mm		Green		
BT-SC2-1636ND7	16 Fr/5.3mm		Orange		
BT-SC2-1836ND7	18 Fr/6.0mm		Red		
BT-SC5-1036ND7	10 Fr/3.3mm	360	Black	Double Swivel Tracheostomy One-way Needleless Irrigation Port MDI Port	10sets/case
BT-SC5-1236ND7	12 Fr/4.0mm		White		
BT-SC5-1436ND7	14 Fr/4.6mm		Green		
BT-SC5-1636ND7	16 Fr/5.3mm		Orange		
BT-SC5-1836ND7	18 Fr/6.0mm		Red		
BT-SC7-1036ND7	10 Fr/3.3mm	360	Black	Double Swivel Tracheostomy Luer Lock Irrigation Port	10sets/case
BT-SC7-1236ND7	12 Fr/4.0mm		White		
BT-SC7-1436ND7	14 Fr/4.6mm		Green		
BT-SC7-1636ND7	16 Fr/5.3mm		Orange		
BT-SC7-1836ND7	18 Fr/6.0mm		Red		
BT-SC8-1036ND7	10 Fr/3.3mm	360	Black	Double Swivel Tracheostomy Luer Lock Irrigation Port MDI Port	10sets/case
BT-SC8-1236ND7	12 Fr/4.0mm		White		
BT-SC8-1436ND7	14 Fr/4.6mm		Green		
BT-SC8-1636ND7	16 Fr/5.3mm		Orange		
BT-SC8-1836ND7	18 Fr/6.0mm		Red		



ACCESSORY CODE:(If Flex tube or Disconnecting Wedge is desired, then include a W, F, WF in Model.#)
BT-SCX-XXXXNWF7:72Hr Stop Valve type+W(Disconnecting Wedge)+F(Flex Tube)