


# "NONSTONE" DOUBLE PIGTAIL URETERAL STENT SET




The Nonstone Ureteral Stent is designed to provide the easy, rapid and comfortable placement benefits of a firm stent. It is generally used for small stones drainage after patients have received an Extracorporeal Shock Wave Lithotripsy operation. Additionally, the stent can function to prevent ureter restenosis after surgery.

- COMFORTABLE APPLICATION*
- SMOOTH MANIPULATION*
- FLEXIBLE RESISTANCE*
- PRECISION POSITIONING*
- QUALITY DESIGN*

## Comfortably Fits Every Anatomy And Every Stone Removal Procedure


- A quality ureteral stent must meet both the procedural needs of the surgeon and accommodate varying preferences.
- **BIOTEQ**'s full line of ureteral stent sets featuring exclusive Thermo Polyurethane construction with high standard of quality and performance to meet the challenges in every different stone removal procedure.
- Selections include Standard Durometer, Enhanced Durometer, Dual Durometer with a selection of curl designs, tip configurations, and surface coatings.

- Unique, inward spiral design to minimize tissue contact for enhanced comfort.
- Dual variable length coil geometry permits balanced stent positioning to minimize the risk of migration.
- Hydrophilic coating provides unequalled surface lubricity which reduces friction to minimize risk of buckling during introduction and placement and to reduce risk of the trauma and encrustation.
- Enhanced Radiopacity for improved visualization.



# STANDARD LOOP STENT SET

**BIOTEQ** Standard Loop Stents are engineered with Enhanced Coil material to provide ease of placement and prevent Stent migration



## Standard Durometer

Flexibility is increased with the standard durometer design to provide maximum comfort needs of the patient.

Used for temporary internal drainage from the ureteropelvic junction to the bladder.

Set includes: Stent, Guidewire, Stent Positioner and Clamp.

Additional option of a suture for repositioning and ease of removal is available.



## Enhanced Durometer

Ease of placement is maximized with the enhanced durometer design without compromising the comfort needs of the patient.

Used for temporary internal drainage from the ureteropelvic junction to the bladder.

Set includes: Stent, Guidewire, Stent Positioner and Clamp.

Additional option of a suture for repositioning and ease of removal is available.

- Enhanced coil material is engineered to provide ease of placement and prevent stent migration
- Thin wall design promotes drainage and patency
- Multiple, large side ports promote drainage
- Attached suture for positioning and subsequent removal without the need for repeat cystoscopy
- Firm, smooth surface eases introduction and radiopacity improves visualization

# STANDARD LOOP URETERAL STENT SET ORDERING INFORMATION

*Standard Durometer/without String	
Item Number	Length (cm)
4.8Fr Pyelo-Close, Bladder-Open	
BT-DJ-4820-110	20
BT-DJ-4822-110	22
BT-DJ-4824-110	24
BT-DJ-4826-110	26
BT-DJ-4828-110	28
BT-DJ-4830-110	30
5Fr Pyelo-Close, Bladder-Open	
BT-DJ-5020-110	20
BT-DJ-5022-110	22
BT-DJ-5024-110	24
BT-DJ-5026-110	26
BT-DJ-5028-110	28
BT-DJ-5030-110	30
6Fr Pyelo-Close, Bladder-Open	
BT-DJ-6020-110	20
BT-DJ-6022-110	22
BT-DJ-6024-110	24
BT-DJ-6026-110	26
BT-DJ-6028-110	28
BT-DJ-6030-110	30
7Fr Pyelo-Close, Bladder-Open	
BT-DJ-7020-110	20
BT-DJ-7022-110	22
BT-DJ-7024-110	24
BT-DJ-7026-110	26
BT-DJ-7028-110	28
BT-DJ-7030-110	30
8Fr Pyelo-Close, Bladder-Open	
BT-DJ-8020-110	20
BT-DJ-8022-110	22
BT-DJ-8024-110	24
BT-DJ-8026-110	26
BT-DJ-8028-110	28
BT-DJ-8030-110	30
4.8Fr Pyelo-Open, Bladder-Open	
BT-DJ-4820-210	20
BT-DJ-4822-210	22
BT-DJ-4824-210	24
BT-DJ-4826-210	26
BT-DJ-4828-210	28
BT-DJ-4830-210	30
5Fr Pyelo-Open, Bladder-Open	
BT-DJ-5020-210	20
BT-DJ-5022-210	22
BT-DJ-5024-210	24
BT-DJ-5026-210	26
BT-DJ-5028-210	28
BT-DJ-5030-210	30
6Fr Pyelo-Open, Bladder-Open	
BT-DJ-6020-210	20
BT-DJ-6022-210	22
BT-DJ-6024-210	24
BT-DJ-6026-210	26
BT-DJ-6028-210	28
BT-DJ-6030-210	30
7Fr Pyelo-Open, Bladder-Open	
BT-DJ-7020-210	20
BT-DJ-7022-210	22
BT-DJ-7024-210	24
BT-DJ-7026-210	26
BT-DJ-7028-210	28
BT-DJ-7030-210	30
8Fr Pyelo-Open, Bladder-Open	
BT-DJ-8020-210	20
BT-DJ-8022-210	22
BT-DJ-8024-210	24
BT-DJ-8026-210	26
BT-DJ-8028-210	28
BT-DJ-8030-210	30

*Standard Durometer/with String	
Item Number	Length (cm)
4.8Fr Pyelo-Close, Bladder-Open	
BT-DJ-4820-111	20
BT-DJ-4822-111	22
BT-DJ-4824-111	24
BT-DJ-4826-111	26
BT-DJ-4828-111	28
BT-DJ-4830-111	30
5Fr Pyelo-Close, Bladder-Open	
BT-DJ-5020-111	20
BT-DJ-5022-111	22
BT-DJ-5024-111	24
BT-DJ-5026-111	26
BT-DJ-5028-111	28
BT-DJ-5030-111	30
6Fr Pyelo-Close, Bladder-Open	
BT-DJ-6020-111	20
BT-DJ-6022-111	22
BT-DJ-6024-111	24
BT-DJ-6026-111	26
BT-DJ-6028-111	28
BT-DJ-6030-111	30
7Fr Pyelo-Close, Bladder-Open	
BT-DJ-7020-111	20
BT-DJ-7022-111	22
BT-DJ-7024-111	24
BT-DJ-7026-111	26
BT-DJ-7028-111	28
BT-DJ-7030-111	30
8Fr Pyelo-Close, Bladder-Open	
BT-DJ-8020-111	20
BT-DJ-8022-111	22
BT-DJ-8024-111	24
BT-DJ-8026-111	26
BT-DJ-8028-111	28
BT-DJ-8030-111	30
4.8Fr Pyelo-Open, Bladder-Open	
BT-DJ-4820-211	20
BT-DJ-4822-211	22
BT-DJ-4824-211	24
BT-DJ-4826-211	26
BT-DJ-4828-211	28
BT-DJ-4830-211	30
5Fr Pyelo-Open, Bladder-Open	
BT-DJ-5020-211	20
BT-DJ-5022-211	22
BT-DJ-5024-211	24
BT-DJ-5026-211	26
BT-DJ-5028-211	28
BT-DJ-5030-211	30
6Fr Pyelo-Open, Bladder-Open	
BT-DJ-6020-211	20
BT-DJ-6022-211	22
BT-DJ-6024-211	24
BT-DJ-6026-211	26
BT-DJ-6028-211	28
BT-DJ-6030-211	30
7Fr Pyelo-Open, Bladder-Open	
BT-DJ-7020-211	20
BT-DJ-7022-211	22
BT-DJ-7024-211	24
BT-DJ-7026-211	26
BT-DJ-7028-211	28
BT-DJ-7030-211	30
8Fr Pyelo-Open, Bladder-Open	
BT-DJ-8020-211	20
BT-DJ-8022-211	22
BT-DJ-8024-211	24
BT-DJ-8026-211	26
BT-DJ-8028-211	28
BT-DJ-8030-211	30

*Enhanced Durometer/without String	
Item Number	Length (cm)
4.8Fr Pyelo-Close, Bladder-Open	
BT-DJ-4820-110-H	20
BT-DJ-4822-110-H	22
BT-DJ-4824-110-H	24
BT-DJ-4826-110-H	26
BT-DJ-4828-110-H	28
BT-DJ-4830-110-H	30
5Fr Pyelo-Close, Bladder-Open	
BT-DJ-5020-110-H	20
BT-DJ-5022-110-H	22
BT-DJ-5024-110-H	24
BT-DJ-5026-110-H	26
BT-DJ-5028-110-H	28
BT-DJ-5030-110-H	30
6Fr Pyelo-Close, Bladder-Open	
BT-DJ-6020-110-H	20
BT-DJ-6022-110-H	22
BT-DJ-6024-110-H	24
BT-DJ-6026-110-H	26
BT-DJ-6028-110-H	28
BT-DJ-6030-110-H	30
7Fr Pyelo-Close, Bladder-Open	
BT-DJ-7020-110-H	20
BT-DJ-7022-110-H	22
BT-DJ-7024-110-H	24
BT-DJ-7026-110-H	26
BT-DJ-7028-110-H	28
BT-DJ-7030-110-H	30
8Fr Pyelo-Close, Bladder-Open	
BT-DJ-8020-110-H	20
BT-DJ-8022-110-H	22
BT-DJ-8024-110-H	24
BT-DJ-8026-110-H	26
BT-DJ-8028-110-H	28
BT-DJ-8030-110-H	30
4.8Fr Pyelo-Open, Bladder-Open	
BT-DJ-4820-210-H	20
BT-DJ-4822-210-H	22
BT-DJ-4824-210-H	24
BT-DJ-4826-210-H	26
BT-DJ-4828-210-H	28
BT-DJ-4830-210-H	30
5Fr Pyelo-Open, Bladder-Open	
BT-DJ-5020-210-H	20
BT-DJ-5022-210-H	22
BT-DJ-5024-210-H	24
BT-DJ-5026-210-H	26
BT-DJ-5028-210-H	28
BT-DJ-5030-210-H	30
6Fr Pyelo-Open, Bladder-Open	
BT-DJ-6020-210-H	20
BT-DJ-6022-210-H	22
BT-DJ-6024-210-H	24
BT-DJ-6026-210-H	26
BT-DJ-6028-210-H	28
BT-DJ-6030-210-H	30
7Fr Pyelo-Open, Bladder-Open	
BT-DJ-7020-210-H	20
BT-DJ-7022-210-H	22
BT-DJ-7024-210-H	24
BT-DJ-7026-210-H	26
BT-DJ-7028-210-H	28
BT-DJ-7030-210-H	30
8Fr Pyelo-Open, Bladder-Open	
BT-DJ-8020-210-H	20
BT-DJ-8022-210-H	22
BT-DJ-8024-210-H	24
BT-DJ-8026-210-H	26
BT-DJ-8028-210-H	28
BT-DJ-8030-210-H	30

*Enhanced Durometer/with String	
Item Number	Length (cm)
4.8Fr Pyelo-Close, Bladder-Open	
BT-DJ-4820-111-H	20
BT-DJ-4822-111-H	22
BT-DJ-4824-111-H	24
BT-DJ-4826-111-H	26
BT-DJ-4828-111-H	28
BT-DJ-4830-111-H	30
5Fr Pyelo-Close, Bladder-Open	
BT-DJ-5020-111-H	20
BT-DJ-5022-111-H	22
BT-DJ-5024-111-H	24
BT-DJ-5026-111-H	26
BT-DJ-5028-111-H	28
BT-DJ-5030-111-H	30
6Fr Pyelo-Close, Bladder-Open	
BT-DJ-6020-111-H	20
BT-DJ-6022-111-H	22
BT-DJ-6024-111-H	24
BT-DJ-6026-111-H	26
BT-DJ-6028-111-H	28
BT-DJ-6030-111-H	30
7Fr Pyelo-Close, Bladder-Open	
BT-DJ-7020-111-H	20
BT-DJ-7022-111-H	22
BT-DJ-7024-111-H	24
BT-DJ-7026-111-H	26
BT-DJ-7028-111-H	28
BT-DJ-7030-111-H	30
8Fr Pyelo-Close, Bladder-Open	
BT-DJ-8020-111-H	20
BT-DJ-8022-111-H	22
BT-DJ-8024-111-H	24
BT-DJ-8026-111-H	26
BT-DJ-8028-111-H	28
BT-DJ-8030-111-H	30
4.8Fr Pyelo-Open, Bladder-Open	
BT-DJ-4820-211-H	20
BT-DJ-4822-211-H	22
BT-DJ-4824-211-H	24
BT-DJ-4826-211-H	26
BT-DJ-4828-211-H	28
BT-DJ-4830-211-H	30
5Fr Pyelo-Open, Bladder-Open	
BT-DJ-5020-211-H	20
BT-DJ-5022-211-H	22
BT-DJ-5024-211-H	24
BT-DJ-5026-211-H	26
BT-DJ-5028-211-H	28
BT-DJ-5030-211-H	30
6Fr Pyelo-Open, Bladder-Open	
BT-DJ-6020-211-H	20
BT-DJ-6022-211-H	22
BT-DJ-6024-211-H	24
BT-DJ-6026-211-H	26
BT-DJ-6028-211-H	28
BT-DJ-6030-211-H	30
7Fr Pyelo-Open, Bladder-Open	
BT-DJ-7020-211-H	20
BT-DJ-7022-211-H	22
BT-DJ-7024-211-H	24
BT-DJ-7026-211-H	26
BT-DJ-7028-211-H	28
BT-DJ-7030-211-H	30
8Fr Pyelo-Open, Bladder-Open	
BT-DJ-8020-211-H	20
BT-DJ-8022-211-H	22
BT-DJ-8024-211-H	24
BT-DJ-8026-211-H	26
BT-DJ-8028-211-H	28
BT-DJ-8030-211-H	30