


BIOTEQ[®] Angiographic Syringe Kit

Procedures made Easier



Bioteque offers custom and OEM capabilities for those seeking complete suite of custom manufacturing services.



We can extrude various tubes to connect with CT/MRI/Angiographic injectors.














Syringes sized to fit a wide variety of injectors including: Liebel-Flarsheim Angiomat[®] 3000 and Angiomat[®] Illumena Injectors...etc.



BIOTEQ[®]

Improving lives

Injector Type	Machine Type	Volume	Syringe Type	Code Number	Diagram	Contents of Kit
CT	Liebel-Flarsheim	200CC	Single Syringe	BT-CT-200-L-I		Angiographic syringe 200ml*1pc+ Fill Tube *1pc+Low Pressure Connector tube (I type)*1pc
			Dual Syringes	BT-CT-200-L-Y		Angiographic syringe 200ml*2pcs+ Fill Tube *2pcs+Low Pressure Connector tube (Y type)*1pc
Angiography	Liebel-Flarsheim	150CC	Single Syringe	BT-A150LC12-Mh120(B)F		Angiographic syringe 150ml+Fill tube x 1pc+ 1200 psi High Pressure Connecting Tube
		150CC	Single Syringe	BT-Angio-150-LF		Angiographic syringe 150ml*1pc+Fill Tube *1pc
Ancillary Devices				CT-1200-Mh120(B)F		Max. 1200 psi Male Luer HP Rotator + 120cm TPU tubing with nylon braided (ID=1.8mm)+Female Luer Lock;
				CT-1200-M120(B)F		Max. 1200 psi, Rotating Male Luer Lock+120cm TPU tubing with nylon braided (ID=1.8mm) +Female Luer Lock;
				CT-1200-Mf120(B)F		Max. 1200 psi Fixed Male Luer Lock + 120cm TPU tubing with nylon braided (ID=1.8mm)+Female Luer Lock;
				BT-150-Y		Coil+Y Type low pressure tube 150cm *1Pc
				BT-150-I		Coil+I Type low pressure tube 150cm *1Pc
				BT-225-Y		Coil+Y Type low pressure tube 225cm *1Pc
				BT-225-Y1		Coil+Y Type low pressure tube 225cm *1Pc+ Check valve