



CENTRAL VENOUS CATHETER SET



When you need a safe and reliable venous access in a person receiving Interventional Radiology procedures, you need the performance of Bioteq Central Venous Catheter Set.

Our products are designed for procedure requiring total precision.

Economical Performance...easier to handle



CVC Features

- Thermo sensitive material allows the tube to soften at body temperature for greater patient comfort.
- Softer catheter tip decrease risk to vascular trauma, radiopaque catheter(BaSO₄) ensures the correct placement of the catheter tip.
- Offer 1cm and 5cm interval depth markings to facilitate catheter placement and clear lumen size and hub colour coding indication.
- Thin wall construction increase flow rates.
- Fixation wing and fastener ensure a complete stability.
- Clearly marked French and gauge size on extension, length and French size on hub.

Central Venous Catheter Set (Single Lumen)

Product No.	Indwelling Catheter	Stainless Steel Guidewire	IV Catheter	Syringe	Dilator	Introducer Needle	Scalpel
BT-CVC-1L-14-08	14Ga x 08cm	0.035" x 45cm	18Ga x 32mm	5cc	8Fr	18Ga x 70mm	#11
BT-CVC-1L-14-13	14Ga x 13cm	0.035" x 45cm	18Ga x 32mm	5cc	8Fr	18Ga x 70mm	#11
BT-CVC-1L-14-16	14Ga x 16cm	0.035" x 45cm	18Ga x 32mm	5cc	8Fr	18Ga x 70mm	#11
BT-CVC-1L-14-20	14Ga x 20cm	0.035" x 68cm	18Ga x 32mm	5cc	8Fr	18Ga x 70mm	#11
BT-CVC-1L-14-30	14Ga x 30cm	0.035" x 68cm	18Ga x 32mm	5cc	8Fr	18Ga x 70mm	#11
BT-CVC-1L-16-08	16Ga x 08cm	0.032" x 45cm	18Ga x 32mm	5cc	7Fr	18Ga x 70mm	#11
BT-CVC-1L-16-13	16Ga x 13cm	0.032" x 45cm	18Ga x 32mm	5cc	7Fr	18Ga x 70mm	#11
BT-CVC-1L-16-16	16Ga x 16cm	0.032" x 45cm	18Ga x 32mm	5cc	7Fr	18Ga x 70mm	#11
BT-CVC-1L-16-20	16Ga x 20cm	0.032" x 68cm	18Ga x 32mm	5cc	7Fr	18Ga x 70mm	#11
BT-CVC-1L-16-30	16Ga x 30cm	0.032" x 68cm	18Ga x 32mm	5cc	7Fr	18Ga x 70mm	#11
BT-CVC-1L-18-08	18Ga x 08cm	0.021" x 45cm	20Ga x 32mm	5cc	6Fr	18Ga x 70mm	#11
BT-CVC-1L-18-13	18Ga x 13cm	0.021" x 45cm	20Ga x 32mm	5cc	6Fr	18Ga x 70mm	#11
BT-CVC-1L-18-16	18Ga x 16cm	0.021" x 45cm	20Ga x 32mm	5cc	6Fr	18Ga x 70mm	#11
BT-CVC-1L-18-20	18Ga x 20cm	0.021" x 68cm	20Ga x 32mm	5cc	6Fr	18Ga x 70mm	#11
BT-CVC-1L-18-30	18Ga x 30cm	0.021" x 68cm	20Ga x 32mm	5cc	6Fr	18Ga x 70mm	#11
BT-CVC-1L-20-08	20Ga x 08cm	0.021" x 45cm	22Ga x 32mm	5cc	5Fr	21Ga x 40mm	#11
BT-CVC-1L-20-13	20Ga x 13cm	0.021" x 45cm	22Ga x 32mm	5cc	5Fr	21Ga x 40mm	#11
BT-CVC-1L-20-16	20Ga x 16cm	0.021" x 45cm	22Ga x 32mm	5cc	5Fr	21Ga x 40mm	#11
BT-CVC-1L-20-20	20Ga x 20cm	0.021" x 68cm	22Ga x 32mm	5cc	5Fr	21Ga x 40mm	#11
BT-CVC-1L-20-30	20Ga x 30cm	0.021" x 68cm	22Ga x 32mm	5cc	5Fr	21Ga x 40mm	#11



Set components:

Central venous catheter x 1, Clamp+Fastener x 1, Spring Wire Guide x 1, Tissue Dilator x1, Scalpel x 1, Raulerson Spring Wire Syringe x 1, Introducer Needle x1, I.V. Catheter x 1

Central Venous Catheter Set (Double Lumen)

Product No.	Indwelling Catheter	Stainless Steel Guidewire	IV Catheter	Syringe	Dilator	Introducer Needle	Scalpel
BT-CVC-2L-70-08	7Fr x 08cm	0.035" x 45cm	18Ga x 32mm	5cc	8.5Fr	18Ga x 70mm	#11
BT-CVC-2L-70-13	7Fr x 13cm	0.035" x 45cm	18Ga x 32mm	5cc	8.5Fr	18Ga x 70mm	#11
BT-CVC-2L-70-16	7Fr x 16cm	0.035" x 45cm	18Ga x 32mm	5cc	8.5Fr	18Ga x 70mm	#11
BT-CVC-2L-70-20	7Fr x 20cm	0.035" x 68cm	18Ga x 32mm	5cc	8.5Fr	18Ga x 70mm	#11
BT-CVC-2L-70-30	7Fr x 30cm	0.035" x 68cm	18Ga x 32mm	5cc	8.5Fr	18Ga x 70mm	#11
BT-CVC-2L-50-08	5Fr x 08cm	0.021" x 45cm	20Ga x 32mm	5cc	6Fr	21Ga x 40mm	#11
BT-CVC-2L-50-13	5Fr x 13cm	0.021" x 45cm	20Ga x 32mm	5cc	6Fr	21Ga x 40mm	#11
BT-CVC-2L-50-16	5Fr x 16cm	0.021" x 45cm	20Ga x 32mm	5cc	6Fr	21Ga x 40mm	#11
BT-CVC-2L-50-20	5Fr x 20cm	0.021" x 68cm	20Ga x 32mm	5cc	6Fr	21Ga x 40mm	#11
BT-CVC-2L-50-30	5Fr x 30cm	0.021" x 68cm	20Ga x 32mm	5cc	6Fr	21Ga x 40mm	#11
BT-CVC-2L-40-08	4Fr x 08cm	0.021" x 45cm	22Ga x 32mm	5cc	5Fr	21Ga x 40mm	#11
BT-CVC-2L-40-13	4Fr x 13cm	0.021" x 45cm	22Ga x 32mm	5cc	5Fr	21Ga x 40mm	#11
BT-CVC-2L-40-16	4Fr x 16cm	0.021" x 45cm	22Ga x 32mm	5cc	5Fr	21Ga x 40mm	#11
BT-CVC-2L-40-20	4Fr x 20cm	0.021" x 68cm	22Ga x 32mm	5cc	5Fr	21Ga x 40mm	#11
BT-CVC-2L-40-30	4Fr x 30cm	0.021" x 68cm	22Ga x 32mm	5cc	5Fr	21Ga x 40mm	#11

Set components:

Central venous catheter x 1, Clamp+Fastener x 1, Spring Wire Guide x 1, Tissue Dilator x1, Scalpel x 1, Raulerson Spring Wire Syringe x 1, Introducer Needle x1, I.V. Catheter x 1



Central Venous Catheter Set (Triple Lumen)

Product No.	Indwelling Catheter	Stainless Steel Guidewire	IV Catheter	Syringe	Dilator	Introducer Needle	Scalpel
BT-CVC-3L-70-08	7Fr x 08cm	0.035" x 45cm	18Ga x 32mm	5cc	8.5Fr	18Ga x 70mm	#11
BT-CVC-3L-70-13	7Fr x 13cm	0.035" x 45cm	18Ga x 32mm	5cc	8.5Fr	18Ga x 70mm	#11
BT-CVC-3L-70-16	7Fr x 16cm	0.035" x 45cm	18Ga x 32mm	5cc	8.5Fr	18Ga x 70mm	#11
BT-CVC-3L-70-20	7Fr x 20cm	0.035" x 68cm	18Ga x 32mm	5cc	8.5Fr	18Ga x 70mm	#11
BT-CVC-3L-70-30	7Fr x 30cm	0.035" x 68cm	18Ga x 32mm	5cc	8.5Fr	18Ga x 70mm	#11
BT-CVC-3L-55-08	5.5Fr x 08cm	0.021" x 45cm	22Ga x 32mm	5cc	6Fr	21Ga x 40mm	#11
BT-CVC-3L-55-13	5.5Fr x 13cm	0.021" x 45cm	22Ga x 32mm	5cc	6Fr	21Ga x 40mm	#11
BT-CVC-3L-55-16	5.5Fr x 16cm	0.021" x 45cm	22Ga x 32mm	5cc	6Fr	21Ga x 40mm	#11
BT-CVC-3L-55-20	5.5Fr x 20cm	0.021" x 68cm	22Ga x 32mm	5cc	6Fr	21Ga x 40mm	#11
BT-CVC-3L-55-30	5.5Fr x 30cm	0.021" x 68cm	22Ga x 32mm	5cc	6Fr	21Ga x 40mm	#11

Set components:

Central venous catheter x 1, Clamp+Fastener x 1, Spring Wire Guide x 1, Tissue Dilator x1, Scalpel x 1, Raulerson Spring Wire Syringe x 1, Introducer Needle x1, I.V. Catheter x 1

