

Guidewire-PTFE Pre-Coated Coil

Each case is different

BIOTEQ® offers **PTFE Pre-Coated Coil Guidewire** in an array of sizes and specifications to suit multiple vascular procedures.

When it comes to safety, everybody counts.

- Pre-coated Coil Guidewire.
- Uniform coating of PTFE provides lubricity and a smooth outer surface for easy insertion.
- Diameters: 0.025", 0.035", 0.038";Lengths: 80cm, 150cm, 260cm.
- Straight and J type available.

Straight

BT-GW-S- α - β ; S stands for tip shape of guidewire is straight.

Guidewire Diameter(α)		Guidewire Length(β)	
α = 018, 021, 025, 035, 038	0.018 inch, 0.021 inch, 0.025 inch, 0.035 inch, 0.038 inch	β = 40, 60, 80, 120, 150, 180, 260, 300	40cm, 60cm, 80cm, 120cm, 150cm, 180cm, 260cm,300cm

J Type

BT-GW-J3- α - β ; J3 stands for tip shape of guidewire is 3.0mm-J shape.

Guidewire Diameter(α)		Guidewire Length(β)	
α = 018, 021, 025, 035, 038	0.018 inch, 0.021 inch, 0.025 inch, 0.035 inch,0.038 inch	β = 40, 60, 80, 120, 150, 260, 300	40cm, 60cm, 80cm, 120cm, 150cm, 260cm,300cm

