

# High Pressure Connecting Tube- High Pressure Lines

## Each case is different

**BIOTEQ® High Pressure Lines** offers a variety of different configurations for contrast injection or pressure transduction and monitoring.

## When it comes to safety, everybody counts

- Made of the medical grade material with pressure resistance, doesn't leak and crack when injection with high-pressure.
- Material includes: Nylon braided TPU, TPU, PVC.
- Intended to access the vasculature for the controlled infusion of diagnostic materials such as contrast media or for the pressure conduction and the monitor uses.
- 1,200, 1,000, 900 psi can be used in a broad range of interventional cardiology applications.
- Overall length: 15~200cm.



## 1200psi

Product Code	Materials	Specifications
CT-1200-Mh120(B)F	TPU tubing with nylon braided	ID:1.8mm, Length 120cm, with Male Luer HP Rotator, Female Luer Lock
CT-1200-Mf100(B)F	TPU tubing with nylon braided	ID:1.8mm, Length 100cm, with Fixed Male Luer Lock, Female Luer Lock
CT-1200-M120(B)F	TPU tubing with nylon braided	ID:1.8mm, Length 120cm, with Rotating Male Luer Lock, Female Luer Lock
CT-1200-M120(T)F	TPU tubing	ID:1.8mm, Length 120cm, with Rotating Male Luer Lock, Female Luer Lock
CT-1200-M120(V1)F	PVC	ID:1.8mm, Length 120cm, with Rotating Male Luer Lock, Female Luer Lock

## 1000psi

Product Code	Materials	Specifications
CT-1000-Mh100(T)F	TPU tubing	ID:1.8mm, Length 100cm, with Male Luer HP Rotator, Female Luer Lock
CT-1000-M100(T)F	TPU tubing	ID:1.8mm, Length 100cm, with Rotating Male Luer Lock, Female Luer Lock

## 900psi

Product Code	Materials	Specifications
CT-900-Mh100F	PVC	ID:2.2mm, Length 100cm, with Male Luer HP Rotator, Female Luer Lock
CT-900-M120(T)F	TPU tubing	ID:2.1mm, Length 120cm, with Rotating Male Luer Lock, Female Luer Lock