

Manifold-2 Ports & 3 Ports

Each case is different

BIOTEQ® Manifold is used to control the flow of the solution and allocate the fluid as a disposable device in PCI procedure.

When it comes to safety, everybody counts.

- 250 psi and 500psi; 2~4 ports ;blue and white handles are available.
- Rotating male luer lock can be connected with other accessories more flexible;
- Off/ on handles are optional.



500psi, Blue and off handles, 3 ports

500psi, Blue and off handles, 3 Ports

Product Code	Specifications
BT-MF-25BF2B-TB	250psi, Blue Handles, Off Handles, 2 Ports, Block Body, Sterile
BT-MF-25WO2B-TB	250psi, White Handles, On Handles, 2 Ports, Block Body, Sterile
BT-MF-50BF3B-TB	500psi, Blue Handles, Off Handles, 3 Ports, Block Body, Sterile