



PTCA Guidewire

Each case is different

BIOTEQ® PTCA Guidewire is intended to facilitate the placement of balloon dilatation catheters or stenting during percutaneous transluminal coronary angioplasty (PTCA), not suitable for coronary chronic total occlusions (CTOs) lesions.

When it comes to safety, everybody counts.

- Coil surface: Hydrophilic coated.
- Dimension: 0.014".
- Length: 195, 300cm.
- Radiopaque tip length: 3cm.
- Tip design: Straight and J.
- Tip stiffness: Standard, Intermediate, Floppy, Extra Floppy.
- Supplied sterile in peel-open packages. Intended for one-time use.



PTCA Guidewire

BT-GW-series-PCI-series	Coil Surface Hydrophilic Coated	Stent Zone Stiffness	Tip Stiffness
BT-GW-S-014-180-PCI	0.014"x195cm, Straight-tip	Standard	Standard
BT-GW-S-014-180-PCI-IN			Intermediate
BT-GW-S-014-180-PCI-FL			Floppy
BT-GW-S-014-180-PCI-EF			Extra Floppy
BT-GW-S-014-180-PCI-D	0.014"x195cm, Straight-tip	Stiff	Standard
BT-GW-S-014-180-PCI-IN-D			Intermediate
BT-GW-S-014-180-PCI-FL-D			Floppy
BT-GW-S-014-180-PCI-EF-D			Extra Floppy
BT-GW-J-014-180-PCI	0.014"x195cm, J-tip	Standard	Standard
BT-GW-J-014-180-PCI-IN			Intermediate
BT-GW-J-014-180-PCI-FL			Floppy
BT-GW-J-014-180-PCI-EF			Extra Floppy
BT-GW-J-014-180-PCI-D	0.014"x195cm, J-tip	Stiff	Standard
BT-GW-J-014-180-PCI-IN-D			Intermediate
BT-GW-J-014-180-PCI-FL-D			Floppy
BT-GW-J-014-180-PCI-EF-D			Extra Floppy
BT-GW-S-014-300-PCI	0.014"x300cm, Straight-tip	Standard	Standard
BT-GW-S-014-300-PCI-IN			Intermediate
BT-GW-S-014-300-PCI-FL			Floppy
BT-GW-S-014-300-PCI-EF			Extra Floppy
BT-GW-S-014-300-PCI-D	0.014"x300cm, Straight-tip	Stiff	Standard
BT-GW-S-014-300-PCI-IN-D			Intermediate
BT-GW-S-014-300-PCI-FL-D			Floppy
BT-GW-S-014-300-PCI-EF-D			Extra Floppy
BT-GW-J-014-300-PCI	0.014"x300cm, J-tip	Standard	Standard
BT-GW-J-014-300-PCI-IN			Intermediate
BT-GW-J-014-300-PCI-FL			Floppy
BT-GW-J-014-300-PCI-EF			Extra Floppy
BT-GW-J-014-300-PCI-D	0.014"x300cm, J-tip	Stiff	Standard
BT-GW-J-014-300-PCI-IN-D			Intermediate
BT-GW-J-014-300-PCI-FL-D			Floppy
BT-GW-J-014-300-PCI-EF-D			Extra Floppy