

Y Connector Set

Each case is different

BIOTEQ® Y Connector Set provides easier way for the doctor to operate in the PCI procedure.

When it comes to safety, everybody counts.

- Screw Type and Push-Pull Type are available.
- Provides single one hand operated procedure which minimizes backflow and blood loss.
- Kit consists of Hemostasis Valve, sideon Female Luer/ Sideon Tubing and Stopcock, Insertion Tool with Large Hub, Red/Copper Torquer.
- Large bore of torque device suits for all sizes of angiographic catheters.



Push-Pull Y Connector, Torquer, Insertion Tool



Screw Type Y Connector, Torquer, Insertion Tool



Screw Type Y Connector, side-on Tubing and Stopcock Torquer, Insertion Tool

Y Connector Set– Push-Pull Type

Product Code	Specifications
BT-HV-7PFLCR-TB	7F, Push-Pull, Sideon Female Luer, Insertion Tool with Large Hub, Red/Copper Torquer
BT-HV-7PPLCR-TB	7F, Push-Pull, Sideon Tubing and Stopcock, Insertion Tool with Large Hub, Red/Copper Torquer
BT-HV-7PFLCW-TB	7F, Push-Pull, Sideon Female Luer, Insertion Tool with Large Hub, White/Copper Torquer
BT-HV-7PPLCW-TB	7F, Push-Pull, Sideon Tubing and Stopcock, Insertion Tool with Large Hub, White/Copper Torquer
BT-HV-7PFLPW-TB	7F, Push-Pull, Sideon Female Luer, Insertion Tool with Large Hub, White/Plastic Torquer
BT-HV-7PPLPW-TB	7F, Push-Pull, Sideon Tubing and Stopcock, Insertion Tool with Large Hub, White/Plastic Torquer

Y Connector Set– Screw Type

Product Code	Specifications
BT-HV-7SFLCR-TB	7F, Screw, Sideon Female Luer, Insertion Tool with Large Hub, Red/Copper Torquer
BT-HV-7SPLCR-TB	7F, Screw, Sideon Tubing and Stopcock, Insertion Tool with Large Hub, Red/Copper Torquer
BT-HV-7SFLCW-TB	7F, Screw, Sideon Female Luer, Insertion Tool with Large Hub, White/Copper Torquer
BT-HV-7SPLCW-TB	7F, Screw, Sideon Tubing and Stopcock, Insertion Tool with Large Hub, White/Copper Torquer
BT-HV-7SFLPW-TB	7F, Screw, Sideon Female Luer, Insertion Tool with Large Hub, White/Plastic Torquer
BT-HV-7SPLPW-TB	7F, Screw, Sideon Tubing and Stopcock, Insertion Tool with Large Hub, White/Plastic Torquer